




Orthopromed Inc.

External Fixator Systems

www.orthopromed.com




CONTENTS

Products	Page
Single-Sided Integrated External Fixator System	3
Annular External Fixator System	13
Combinational External Fixator System	26
Common External Fixator System	51
Accessories Of External Fixator System	61

>> SINGLE-SIDED INTEGRATED EXTERNAL FIXATOR SYSTEM


Mini Fragment External Fixator-Type A
 REF: 6100-0101
 Matching bone screws: Ø2x70mm 4 pcs
 Matching instruments: 3mm hex wrench, 2mm screwdriver
 Length: 100mm




Mini Fragment External Fixator-Type B
 REF: 6100-0102
 Matching bone screws: Ø2x70mm 4 pcs
 Matching instruments: 3mm hex wrench, 2mm screwdriver
 Length: 90mm



Mini Fragment External Fixator-Type C
 REF: 6100-0103
 Matching bone screws: Ø2.5x50mm 4 pcs
 Matching instruments: 3mm hex wrench, 2mm screwdriver
 Length: 74mm



Mini Fragment External Fixator-Type D
 REF: 6100-0104
 Matching bone screws: Ø2x50mm 4 pcs
 Matching instruments: 3mm hex wrench, 2mm screwdriver
 Length: 45mm



» SINGLE-SIDED INTEGRATED EXTERNAL FIXATOR SYSTEM


Mini Fragment External Fixator-Type E

REF: 6100-0105
Matching bone screws: Ø2x70mm 4 pcs
Matching instruments: 3mm hex wrench, 2mm screwdriver
Length: 60mm




Mini Fragment External Fixator-Type F

REF: 6100-0106
Matching bone screws: Ø2x70mm 4 pcs
Matching instruments: 3mm hex wrench, 2mm screwdriver
Length: 70mm



Mini Fragment External Fixator-Type G

REF: 6100-0107
Matching bone screws: Ø2x70mm 4 pcs
Matching instruments: 3mm hex wrench, 2mm screwdriver
Length: 90mm



Mini Fragment External Fixator-Type H


REF: 6100-0108
Matching bone screws: Ø2x70mm 4 pcs
Matching instruments: 3mm hex wrench, 2mm screwdriver
Length: 70mm



» SINGLE-SIDED INTEGRATED EXTERNAL FIXATOR SYSTEM


Wrist External Fixator-Type A

REF: 6100-0201
Matching bone screws: Ø3.5x80mm 4 pcs
Matching instruments: 3.5mm screwdriver, 4mm hex wrench




Wrist External Fixator-Type B

REF: 6100-0202
Matching bone screws: Ø2.5x60mm 2 pcs, Ø3.5x80mm 2 pcs
Matching instruments: 10mm open spanner, 2.5&3.5mm screwdriver, 4mm hex wrench



Wrist External Fixator-Type D

REF: 6100-0204
Matching bone screws: Ø2.5x60mm 2 pcs, Ø3.5x80mm 2 pcs
Matching instruments: 3.5mm screwdriver, 4mm hex wrench



» SINGLE-SIDED INTEGRATED EXTERNAL FIXATOR SYSTEM


Wrist External Fixator-Type E

REF: 6100-0205
Matching bone screws: Ø2.5x60mm 2 pcs, Ø3.5x80mm 2 pcs
Matching instruments: 3.5mm screwdriver, 4mm hex wrench




Wrist External Fixator-Type F

REF: 6100-0206
Matching bone screws: Ø3.5x80mm 4 pcs,
Matching instruments: 3.5mm screwdriver, 5mm hex wrench,
10mm open spanner




Wrist External Fixator (Radiolucent Type)

REF: 6100-0207
Matching bone screws: Ø2.5x60mm 2 pcs, Ø3.5x80mm 2 pcs
Matching instruments: 3.5mm screwdriver, 4mm hex wrench



Wrist External Fixator (Radiolucent Type)

REF: 6100-0208
Matching bone screws: Ø2.5x60mm 2 pcs, Ø3.5x80mm 2 pcs
Matching instruments: 3.5mm screwdriver, 5mm hex wrench



» SINGLE-SIDED INTEGRATED EXTERNAL FIXATOR SYSTEM

Dynamic Axial External Fixator(Straight Type)

XL
REF: 6100-0301
Matching bone screws: Ø6x150mm 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

L
REF: 6100-0302
Matching bone screws: Ø6x150mm 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

M
REF: 6100-0303
Matching bone screws: Ø6x130mm 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

S
REF: 6100-0304
Matching bone screws: Ø5x120mm 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

XS
REF: 6100-0305
Matching bone screws: Ø4x100mm 4 pcs
Matching instruments: 4mm hex wrench, 4mm screwdriver



Dynamic Axial External Fixator(T Type)

XL
REF: 6100-0401
Matching bone screws: Ø6x150mm 2 pcs, HBØ6.0x150mm 2 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

L
REF: 6100-0402
Matching bone screws: Ø6x150mm 2 pcs, HBØ6.0x150mm 2 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

M
REF: 6100-0403
Matching bone screws: Ø6x130mm 2 pcs, HBØ6.0x150mm 2 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

S
REF: 6100-0404
Matching bone screws: Ø6x120mm 2 pcs, HBØ5.0x130mm 2 pcs
Matching instruments: 5mm hex wrench, 5mm screwdriver

XS
REF: 6100-0405
Matching bone screws: Ø4x100mm 4 pcs
Matching instruments: 4mm hex wrench, 4mm screwdriver



» SINGLE-SIDED INTEGRATED EXTERNAL FIXATOR SYSTEM

Dynamic Axial External Fixator(Ankle Joint Type)

XL
REF: 6100-0501
Matching bone screws: Ø6x150mm 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

L
REF: 6100-0502
Matching bone screws: Ø6x150mm 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

M
REF: 6100-0503
Matching bone screws: Ø6x130mm 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver



Dynamic Axial External Fixator(Hip Joint Type)

XL
REF: 6100-0601
Matching bone screws: Ø6x150mm 2 pcs, HBØ6.0x150mm 2 pc
Matching instruments: 6mm hex wrench, 6mm screwdriver

L
REF: 6100-0602
Matching bone screws: Ø5x150mm 2 pcs, HBØ6.0x150mm 2 pc
Matching instruments: 6mm hex wrench, 6mm screwdriver




Petrochanteric External Fixator

REF: 6100-0701
Matching bone screws: Ø6x150mm 2 pcs, HBØ6.0x180mm 2 pcs
Matching instruments: 5mm hex wrench, 6mm screwdriver

» SINGLE-SIDED INTEGRATED EXTERNAL FIXATOR SYSTEM

Elbow Joint External Fixator

REF: 6100-0801
Matching bone screws: Ø5x110mm 4 pcs
Matching instruments: 5mm hex wrench, 3mm hex wrench, 6mm screwdriver



Knee Joint External Fixator


REF: 6100-0901
Matching bone screws: Ø6x150mm 2 pcs, Ø6x130mm 2pcs
Matching instruments: 6mm hex wrench, 3mm hex wrench, 6mm screwdriver



» SINGLE-SIDED INTEGRATED EXTERNAL FIXATOR SYSTEM

Pelvic External Fixator-Type A

REF: 6100-1001
Matching bone screws: HBØ6x150mm 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver
Option: 5mm Pin




Pelvic External Fixator-Type B

REF: 6100-1002
Matching bone screws: HBØ6x150mm 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver
Option: 5mm Pin



Pelvic External Fixator-Type C

REF: 6100-1003
Matching bone screws: HBØ6x150mm 4 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver
Option: 5mm Pin



» SINGLE-SIDED INTEGRATED EXTERNAL FIXATOR SYSTEM

Bone Lengthening External Fixator-Type A

REF: 6100-1101(250mm) For Children
Matching bone screws: Ø4.0x120mm 6 pcs
Matching instruments: 4mm hex wrench, 4mm screwdriver

REF: 6100-1102(300mm)
Matching bone screws: Ø6x150mm 6 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

REF: 6100-1103(350mm)
Matching bone screws: Ø6x150mm 6 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver

REF: 6100-1104(400mm)
Matching bone screws: Ø6x150mm 6 pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver




Bone Lengthening External Fixator-Type B

REF: 6100-1105(250mm)
Matching bone screws: Ø6x150mm 4 pcs (Select 5mm bone screw)
Matching instruments: 6mm hex wrench, 6mm screwdriver

REF: 6100-1106(300mm)
Matching bone screws: Ø6x150mm 4 pcs (select 5mm bone screw)
Matching instruments: 6mm hex wrench, 6mm screwdriver

Bone Lengthening External Fixator-Type C

REF: 6100-1107(300mm)
Matching bone screws: Ø6x150mm 6 pcs,
HBØ6.0x150mm 2pcs
Matching instruments: 6mm hex wrench, 6mm screwdriver



» SINGLE-SIDED INTEGRATED EXTERNAL FIXATOR SYSTEM

Guide Rail

REF: Z04-01
Size: 200, 250, 300, 400mm



Intermediate Clamp

REF: Z04-02



Outer End Clamp

REF: Z04-03



Adjustment Device

REF: Z04-04
Size: 50, 70, 110mm



T-shaped Correction Clamp

REF: Z04-05
Suitable for correcting foot eversion Or varus



T-clamp

REF: Z04-06



Elastic Push Block

REF: Z04-07
Apply pressure to the fracture line




Double Ring Bridge Block

REF: Z04-08
For bridging with ring frame




Fine Tuning Rotary Clamp

REF: Z04-09
For correction of angular deformity



Multi-directional Clamp

REF: Z04-10
Suitable for correction of deformed rotation and movement



Universal T-clamp


REF: Z04-11
For joint end



» ANNULAR EXTERNAL FIXATOR SYSTEM

Transparent Ring Frame Combination (PEEK)

REF: 6100-1204



Transparent Ring Frame Combination (PEEK)

REF: 6100-1201

Tibial Fixator Transparant (PEEK)

REF: 6100-1202



Calcaneal Fixator Transparant (PEEK)

REF: 6100-1203



» ANNULAR EXTERNAL FIXATOR SYSTEM

Ring Connection
REF: HH-01



1/3 Rings
REF: HH-02
Size: Ø125, Ø155, Ø175,
Ø195, Ø215



2/3 Rings
REF: HH-03
Size: Ø125, Ø155, Ø175,
Ø195, Ø215



Foot Rings
REF: HH-04



Rod To Rod Coupling I
REF: HH-05
Size: Ø6/6



Rod To Rod Coupling II
REF: HH-06
Size: Ø6/6



Pin To Ring Coupling
REF: HH-07
Size: Ø1.5/3



Two-hole K-wire Pin
REF: HH-10



90° Extension Plate
REF: HH-11



Single Screw Clamp
REF: HH-12



Diagonal Tie
REF: HH-08



Single Direction Joint With Displacement
REF: HH-13



» ANNULAR EXTERNAL FIXATOR SYSTEM

Ring Bar
REF: HH-09
Size: 50, 100, 150



K-wire Pin
REF: HH-14



Quick Connecting Rod
REF: HH-16



Single Direction Joint
REF: HH-15



Micro-motion Part For Rings
REF: HH-18



Ring Clamp
REF: HH-17



T Wrench
REF: HGJ01



Kirschner Wire Tensioning Pliers
REF: HGJ02



Kirschner Wire Guider
REF: HGJ03



Open Spanner
REF: HGJ04
Size: #10



» ANNULAR EXTERNAL FIXATOR SYSTEM

Ring External Fixator

REF: 6100-1201

Used for: Tibia, Femur

Configuration:

Half rings 6 pcs, M6x300 threaded rod 3 pcs, M6x20 thread 6 pcs, M6x20 nut 6 pcs, Pin clamp 16 pcs




2/3 Ring External Fixator

REF: 6100-1202

Used for: Tibia, Femur

Configuration:

Rings with 6 groove 3 pcs, M6x400 threaded rod 3 pcs, Pin clamp 12 pcs



Lower Limb Lengthening External Fixator

REF: 6100-1203

Used for: Tibia, Femur

Configuration:

Rings 4 pcs, Inner threaded connecting rod 4 pcs, M6x180 threaded rod 4 pcs, M6x120 threaded rod 4 pcs, Pin clamp 16 pcs

» ANNULAR EXTERNAL FIXATOR SYSTEM


Ring External Fixator

REF: 6100-1204

Used for: Tibia, Femur

Configuration:

Rings 4 pcs, M6x400 threaded rod 4 pcs, adjusting nut 4 pcs, Pin clamp 16 pcs, Annular tubes 4 pcs



Elbow Joint External Fixator

REF: 6100-1205

Used for: Elbow joint fixator

Configuration:

Rings 4 pcs, Annular diagonal tie 1 pc, single joint 2 pcs, M6x200 threaded rod 8 pcs, Pin clamp 16 pcs, 2-hole connecting plate 1 pc

Ankle External Fixator

REF: 6100-1206

Used for: Equinus

Configuration:

Rings 4 pcs, L connecting 2 pcs, half rings 1 pc, Annular diagonal tie 1 pc, inner threaded post 4 pcs, single joint 2 pcs, M6 threaded rod 8 pcs, Pin clamp 24 pcs



» ANNULAR EXTERNAL FIXATOR SYSTEM

Ring External Fixator

REF: 6100-1207

Used for: Tibia, Femur

Configuration:

Rings with 1 ear, rings with six ears, standard post 6 pcs



Knee Joint External Fixator

REF: 6100-1208

Used for: Knee joint

Configuration:

Rings 4 pcs, foot ring 1 pc, half rings 1 pc, Annular diagonal tie 1 pc, post with 4 edges 2 pcs, single joint 2 pcs, Pin clamp 16 pcs, M6 threaded rod 12 pcs



Ankle Joint External Fixator

REF: 6100-1209

Used for: Equinus

Configuration:

Rings 2 pcs, foot ring 1 pc, L connecting 2 pcs, half rings 1 pc, Annular diagonal tie, Single joint 2 pc, M6x200 threaded rod 3 pcs, Pin clamp 14 pcs

» ANNULAR EXTERNAL FIXATOR SYSTEM

Parts & Instruments Of Ring External Fixator

Ring With Six Ears

REF: 6100-1210 L

REF: 6100-1211 M

REF: 6100-1212 S

Material: Aluminum



Ring With 3 Ears

REF: 6100-1213 L

REF: 6100-1214 M

REF: 6100-1215 S

Material: Aluminum



Half Ring

REF: HTLL0301 L

REF: HTLL0302 M

REF: HTLL0303 S

Material: Aluminum



Foot Ring

REF: HTLL04 L

Material: Aluminum



Post With 4 Edges

REF: 6100-1216 1 holes

REF: 6100-1217 2 holes

REF: 6100-1218 3 holes

REF: 6100-1219 4 holes

REF: 6100-1220 5 holes

REF: 6100-1221 6 holes


Material: Stainless steel



Post With 6 Edges

REF: 6100-1317

Material: Stainless steel



Rod To Ring Bolt

REF: 6100-1230

Material: Stainless steel



Bush

REF: 6100-1222 Ø 4

REF: 6100-1223 Ø 5

REF: 6100-1224 Ø 6

Material: Aluminum



Pin To Ring Coupling

REF: 6100-1225 Ø 1.5

REF: 6100-1226 Ø 2

REF: 6100-1227 Ø 3

Material: Stainless steel



Pin To Post Coupling

REF: 6100-1228 Ø 4mm

REF: 6100-1229 Ø 5mm

Material: Stainless steel



Standard Rod

REF: 6100-1231 L

REF: 6100-1232 M

REF: 6100-1233 S



Packaged Rod

REF: 6100-1234 L

REF: 6100-1235 M

REF: 6100-1236 S



» ANNULAR EXTERNAL FIXATOR SYSTEM

Parts & Instruments Of Ring External Fixator


Rings
REF: 6100-1239 Ø100mm
REF: 6100-1240 Ø120mm
REF: 6100-1241 Ø140mm
REF: 6100-1242 Ø160mm
REF: 6100-1243 Ø180mm
REF: 6100-1244 Ø200mm
REF: 6100-1245 Ø220mm
Material: Aluminum




Half Rings
REF: 6100-1247 Ø100mm
REF: 6100-1248 Ø120mm
REF: 6100-1249 Ø140mm
REF: 6100-1250 Ø160mm
REF: 6100-1251 Ø180mm
REF: 6100-1252 Ø200mm
REF: 6100-1253 Ø220mm
Material: Aluminum



2/3 Rings
REF: 6100-1255 Ø100mm
REF: 6100-1256 Ø120mm
REF: 6100-1257 Ø140mm
REF: 6100-1258 Ø160mm
REF: 6100-1259 Ø180mm
REF: 6100-1260 Ø200mm
REF: 6100-1261 Ø220mm
Material: Aluminum



Foot Rings
REF: 6100-1262 L
REF: 6100-1263 M
REF: 6100-1264 S
Material: Stainless steel



Annular Diagonal Tie
REF: 6100-1237



Threaded Rod M6
REF: 1001-0601/0613(40-400mm)



Adjusting Nut
REF: 6100-1265



Annular Tubes
REF: 6100-1266



Spring
REF: 6100-1267



Single Joint
REF: 6100-1268
Material: Stainless steel



Single Joint
REF: 6100-1269
Material: Stainless steel



Universal Joint
REF: 6100-1270
Material: Stainless steel



» ANNULAR EXTERNAL FIXATOR SYSTEM

Parts & Instruments Of Ring External Fixator

Threaded Column
REF: 6100-1271 1 holes
REF: 6100-1272 2 holes
REF: 6100-1273 3 holes
REF: 6100-1274 4 holes
REF: 6100-1275 5 holes
REF: 6100-1276 6 holes
REF: 6100-1277 7 holes
REF: 6100-1278 8 holes
Material: Stainless steel



Nut Column
REF: 6100-1279 1 holes
REF: 6100-1280 2 holes
REF: 6100-1281 3 holes
REF: 6100-1282 4 holes
REF: 6100-1283 5 holes
REF: 6100-1284 6 holes
REF: 6100-1285 7 holes
REF: 6100-1286 8 holes
Material: Stainless steel



L Connecting Plate
REF: 6100-1287 5 holes
REF: 6100-1288 6 holes
REF: 6100-1289 7 holes
REF: 6100-1290 8 holes
REF: 6100-1291 9 holes
REF: 6100-1292 10 holes
REF: 6100-1293 11 holes
REF: 6100-1294 12 holes
Material: Stainless steel



Connecting Plate
REF: 6100-1295 2 holes
REF: 6100-1296 3 holes
REF: 6100-1297 4 holes
REF: 6100-1298 5 holes
REF: 6100-1299 6 holes
REF: 6100-1300 7 holes
REF: 6100-1301 8 holes
Material: Stainless steel



CrossConnecting Plate
REF: 6100-1302 2 holes
REF: 6100-1303 3 holes
REF: 6100-1304 4 holes
REF: 6100-1305 5 holes
REF: 6100-1306 6 holes
REF: 6100-1307 7 holes
REF: 6100-1308 8 holes
Material: Stainless steel



Inner Threaded Connecting Rod
REF: 6100-1309 80mm
Number: 6100-1310 120mm
Material: Stainless steel



M6 Bolt
REF: HYLL17
Size: M6x10, 12, 16, 20, 25, 30
Material: Stainless steel



Pin Clamp
REF: 6100-1311 Ø4
6100-1312 Ø5-6
Material: Stainless steel



Pin Clamp Thick
REF: 6100-1313 Ø1.5mm
REF: 6100-1314 Ø2.5mm
Material: Stainless steel




Drill Bit Stop Unit
REF: 6100-1315
Material: Stainless steel



Hex Nut M6
REF: 6100-1316
Material: Stainless steel



Wrist Joint External Fixator
REF: HOT001
Used for: wrist
Configuration:
2pcs Ø160 full ring, 3pcs Ø160 half ring, 2pcs hinge, 1pc retractor,
7pcs M6x250 threaded rod




Ankle Joint External Fixator
REF: HOT002
Used for: Ankle joint
Configuration:
4pcs extender, 6pcs M6x150 threaded rod, 3pcs porous plate with 1 hole,
3pcs universal joint, 2pcs Ø160 ring, 2pcs Ø160U ring



Parts & Instruments

Ring
REF: HOTL01
Size:
Ø120, Ø140, Ø150, Ø160,
Ø180, Ø200, Ø220, Ø240
Material: Aluminum



Half Ring
REF: HOTL02
Size:
Ø120, Ø140, Ø150, Ø160, Ø180,
Ø200, Ø220, Ø240
Material: Aluminum




U Ring
REF: HOTL03
Size:
Ø120, Ø140, Ø150, Ø160
Ø180, Ø200
Material: Aluminum




90° Ring
REF: HOTL04
Size:
R90 long, R90 Short, R120 Long, R120 Short
Material: Aluminum




Porous Plate
REF: HOTL05
Size:l, 2, 3, 4, 5, 6, 7, 8 Holes
Material: Stainless steel




Six Corner Connection
REF: HOTL06
Size: 15mm, 30mm
Material: Stainless steel



Connecting Plate
REF: HOTL07
Size: 20, 30, 40, 50
Material: Stainless steel



U Type Ring Extension Rod
REF: HOTL08
Size: 3, 4, 5, 6, 7, 8 Holes
Material: Aluminum



Ball Head Screw
REF: HOTL09
Size: Ø100, Ø150, Ø200, Ø280
Material: Stainless steel



90° Conversion
REF: HOTL10
Material: Stainless steel



Square Nuts
REF: HOTL11
Material: Stainless steel



Nylon Square Nut
REF: HOTL12
Material: Nylon




Cone Washer
REF: HOTL13
Material: Stainless steel



Cone Washer
REF: HOTL14
Material: Stainless steel






Annular Tubes
REF: HOTL15
Size: 10, 20, 30, 40, 50




» ANNULAR EXTERNAL FIXATOR SYSTEM


Parts & Instruments

Connecting Rod REF: HOTL16 Material: Aluminum 	Quick Coupling Rod REF: HOTL17 Material: Stainless steel 	Universal Joint REF: HOTG03 Material: Stainless steel 
Retractor REF: HOTL18 Size: 100, 150, 200 Material: Stainless steel 	Retractor II REF: HOTL19 Material: Stainless steel 	Threaded Rod REF: HYLL05 
Needle Clamp REF: HOTJ01 Material: Stainless steel 	Square Needle Clamp REF: HOTJ02 Material: Stainless steel 	Pin To Ring Coupling REF: HOTJ03 Material: Stainless steel 
Hinge REF: HOTG01 Material: Stainless steel 	L Type Hinge REF: HOTG02 Material: Stainless steel 	Extender REF: HOTY01 Size: 70, 100, 150, 200 Material: Aluminum 
Tight Wire Clamp REF: HOTGJ01 	Socket Wrench #10 REF: HOTGJ02 	Straight Screwdriver #3 REF: HOTGJ03 
Screw Driver REF: HOTGJ04 Size: 6mm 	Double Ended Spanner REF: HOTGJ05 Size: #10 	Socket Wrench REF: HOTGJ06 Size: #10 

» ANNULAR EXTERNAL FIXATOR SYSTEM




Lateral Limb Reconstruction External Fixator
REF: B01
Size: Ø4.0x100, Ø4.0x160



Large Skin Retractor
REF: Q01
Size: Ø2.0mmx250
Ø2.5mmx250



Dual-purpose Skin Retractor
REF: Q03



Small Skin Retractor
REF: Q02
Size: Ø1.5mmx250
Ø2.0mmx250

COMBINATIONAL EXTERNAL FIXATOR SYSTEM

Small Fragment External Fixator Instrument Set 7100-01

PRODUCT NO.	NO.	PRODUCT NAME	SPECIFICATION	QTY
7100-0101	1	Rod To Rod Coupling	Ø5/5	12
7100-0102	2	PinTo Rod Coupling	Ø5/3-4	12
7100-0103	3	4-hole Pin Clamp	Ø5/3-4	3
7100-0104	4	Peri-articular Pin Clamp	Ø3/4	1
7100-0105	5	Post	Ø5/0°	2
7100-0106	6	Post	Ø5/30°	4
7100-0107	7	Apex Pin Self-drilling	HA 3.0x80	4
7100-0108			HA 3.0x100	4
7100-0109			HA 4.0x80	1
7100-0110			HA 4.0x100	4
7100-0111			HA 4.0x120	4
7100-0112	8	Drill Sleeves	Ø3	1
7100-0113			Ø4	1
7100-0114	9	Carbon Fiber Connecting Rod	Ø5x120	4
7100-0115			Ø5x150	4
7100-0116			Ø5x180	2
7100-0117			Ø5x200	4
7100-0118	10	T Wrench	#5	1
7100-0119	11	Stabilization/Reduction Wrench	#5	1
7100-0120	12	Elbow Joint Movement Device	Ø5	1
7100-0121	13	Apex Pin Chuck	Ø3	1
7100-0122			Ø4	1
7100-0123	14	T-shaped Pin Wrench	Ø3	1
7100-0124			Ø4	1
7100-0125	15	Hand Drill		1
7100-0126	16	Thumb Wheel	#5/7	1
7100-0127	17	Aluminum Box		1

COMBINATIONAL EXTERNAL FIXATOR SYSTEM

Part & Instruments Of Small Fragment External Instrument Set

Pin To Rod Coupling REF: 7100-0102 Size: Ø5/3-4 Material: Aluminum 	Rod To Rod Coupling REF: 7100-0101 Size: Ø5/5 Material: Aluminum 	Pin To Rod Coupling REF: ZSSP01 Size: Ø5/3-4 Material: PEEK 
Rod To Rod Coupling REF: ZSSR01 Size: Ø5/5 Material: PEEK 	Pin To Rod Integrated Coupling REF: ZSSP04 Size: Ø3/4/5 Material: Aluminum 	Peri-articular Pin Clamp REF: 7100-0104 Size: Ø5/3-4 Material: Aluminum 
4-hole Pin Clamp REF: 7100-0103 Size: Ø5/3-4 Material: Aluminum 	Pin To Rod Coupling REF: 7100-0127 Size: Ø3/2 Material: Aluminum 	Rod To Rod Coupling REF: 7100-0128 Size: Ø3/3 Material: Aluminum 

>> COMBINATIONAL EXTERNAL FIXATOR SYSTEM

Part & Instruments Of Small Fragment External Instrument Set

Tube To Rod Coupling
REF: 7100-0129
Size: Ø15/5
Material: Aluminum



Connecting Rod
REF: 7100-0130/7100-0114
Material: Aluminum/Carbon fiber




Straight Post
REF: 7100-0105
Size: Ø5/0°
Material: Stainless steel



Elbow Joint Movement Device
REF: 7100-0120
Size: Ø5
Material: Stainless steel



30° post
REF: 7100-0106
Size: Ø5/30°
Material: Stainless steel



Wrist Joint Movement Device
REF: 7100-0601
Size: Ø5
Material: Stainless steel



Compression/distraction Tube
REF: 7100-1101
Size: 15mm
Material: Aluminum



>> COMBINATIONAL EXTERNAL FIXATOR SYSTEM

Part & Instruments Of Small Fragment External Instrument Set

Stabilization/reduction Wrench
REF: 7100-0119
Size: #15




Hand Drill
REF: 7100-0125




Wrench
REF: 7100-0133
Size: #3/#4/#5




T Wrench
REF: 7100-0118
Size: #5



T-shaped Pin Wrench
REF: 7100-0123/7100-0124
Size: Ø3/Ø4



Drill Sleeves
REF: 7100-0112
Size: Ø3, Ø4



Apex Pin Chuck
REF: 7100-0121/7100-0122
Size: Ø3/Ø4



Thumb Wheel
REF: 7100-0126
Size: #5/7




» COMBINATIONAL EXTERNAL FIXATOR SYSTEM

Ulna & Radius External Fixator

REF: 7100-03
Used for: Ulna, Radius


REF	PRODUCT NAME	SPECIFICATION	QTY
7100-0301	Rod To Rod Coupling	Ø5/5	4
7100-0302	4-hole Pin Clamp	Ø5/3-4	2
7100-0303	Post	Ø5/30°	4
7100-0304	Connecting rod	Ø5x200mm	2
7100-0305	Apex Pin Self-drilling	Ø3x120mm	4



Forearm External Fixator

REF: 7100-04
Used for: Forearm


REF	PRODUCT NAME	SPECIFICATION	QTY
7100-0401	Rod To Rod Coupling	Ø5/5	2
7100-0402	4-hole Pin Clamp	Ø5/3-4	2
7100-0403	Post	0°	2
7100-0404	Connecting rod	Ø5x200mm	1
7100-0405	Apex Pin Self-drilling	Ø3x120mm	4



Distal Radius External Fixator

REF: 7100-05
Used for: Ulna, Radius

REF	PRODUCT NAME	SPECIFICATION	QTY
7100-0501	Rod To Rod Coupling	Ø5/5	1
7100-0502	Pin To Rod Clamp	Ø5/3-4	2
7100-0503	Peri-Articular Pin Clamp	Ø5/3-4	1
7100-0504	Connecting rod	Ø5x200mm	1
7100-0505	Apex Pin Self-drilling	Ø3x120mm	4




» COMBINATIONAL EXTERNAL FIXATOR SYSTEM

Wrist Joint External Fixator

REF: 7100-06
Used for: Wrist Joint


REF	PRODUCT NAME	SPECIFICATION	QTY
7100-0601	Wrist Joint Movement Device	Ø5	1
7100-0602	Pin To Rod Coupling	Ø5/3-4	4
7100-0603	Apex Pin Self-drilling	Ø3x120mm	4



Elbow Joint External Fixator

REF: 7100-07
Used for: Elbow Joint


REF	PRODUCT NAME	SPECIFICATION	QTY
7100-0701	Elbow Joint Movement Device	Ø5	1
7100-0702	Pin To Rod Coupling	Ø5/3-4	4
7100-0703	Apex Pin Self-drilling	4x120mm	4



Finder & Palm External Fixator

REF: 7100-08
Used for: Finger, Palm

REF	PRODUCT NAME	SPECIFICATION	QTY
7100-0801	Pin To Rod Coupling	Ø3/2	2
7100-0802	Connecting Rod	Ø3x80mm	1
7100-0803	Apex Pin Self-drilling	Ø2x70mm	4



» COMBINATIONAL EXTERNAL FIXATOR SYSTEM

Humeral Shaft External Fixator

REF: 7100-09
Used for: Humeral Shaft


REF	PRODUCT NAME	SPECIFICATION	QTY
7100-0901	Rod To Rod Coupling	Ø5/5	2
7100-0902	4-hole Pin Clamp	Ø5/3-4	2
7100-0903	Post	Ø5/30°	2
7100-0904	Connecting Rod	Ø5x200mm	1
7100-0905	Apex Pin Self-drilling	Ø4x120mm	4



Proximal Humerus External Fixator

REF: 7100-10
Used for: Proximal Humerus, Condyle Neck

REF	PRODUCT NAME	SPECIFICATION	QTY
7100-1001	Rod To Rod Coupling	Ø5/5	4
7100-1002	Pin To Rod Coupling	Ø5/3-4	2
7100-1003	4-hole Pin Clamp	Ø5/3-4	1
7100-1004	Connecting Rod	Ø5x250mm	2
7100-1005	Connecting Rod	Ø5x200mm	1
7100-1006	Post	Ø5/30°	2
7100-1007	Apex Pin Self-drilling	Ø4x120mm	4




Humerus External Fixator

REF: 7100-11
Used for: Bone Shaft, Distal & Proximal Bone

REF	PRODUCT NAME	SPECIFICATION	QTY
7100-1101	Compression/Distracton Tube	Ø15	1
7100-1102	Tube To Rod Coupling	Ø15/5	2
7100-1103	4-hole Pin Clamp	Ø5/3-4	2
7100-1104	Post	Ø5/30°	2
7100-1105	Apex Pin Self-drilling	Ø4x120mm	4

» COMBINATIONAL EXTERNAL FIXATOR SYSTEM

Large Frament External Fixator Instrument Set 7100-02



PRODUCT NO.	NO.	PRODUCT NAME	SPECIFICATION	QTY
7100-0201	1	Rod To Rod Coupling	Ø8/8	12
7100-0202	2	Pin To Rod Coupling	Ø8/5-6	12
7100-0203	3	5-hole Pin Clamp	Ø8/5-6	3
7100-0204	4	Post	Ø8/0°	2
7100-0205	5	Post	Ø5/30°	4
7100-0206	6	Apex Pin Self-drilling	HA 5.0x120	4
7100-0207			HA 5.0x150	8
7100-0208			HA 5.0x180	4
7100-0209	7	Drill Sleeves	Ø5	2
7100-0210	8	Carbon Fiber Connecting Rod	Ø8x200	4
7100-0211			Ø8x250	4
7100-0212			Ø8x300	4
7100-0213	9	Semi-Circular Curved Rod	Ø8 180	2
7100-0214			Ø8/200	2
7100-0215	10	Stabilization/Reduction Wrench	#7	1
7100-0216	11	Torque Wrench	Ø7	1
7100-0217	12	Thumb Wheel	#5/7	1
7100-0218	13	Apex Pin Chuck	Ø5	1
7100-0219	14	Hand Drill		1
7100-0220	15	T-shaped Pin Wrench	Ø5	1
7100-0221	16	T Wrench	#7	1
7100-0222	17	Aluminum Box		1

COMBINATIONAL EXTERNAL FIXATOR SYSTEM

Parts & Instruments Of Large Fragment External Fixator Instrument Set

Pin To Rod Coupling REF: 7100-0201 Size: Ø8/5-6 Material: Aluminum 	Rod To Rod Coupling REF: 7100-0202 Size: Ø8/8 Material: Aluminum 	Pin To Rod Coupling REF: ZSSP01 Size: Ø8/5-6 Material: Aluminum 
Rod To Rod Coupling REF: ZSSR01 Size: Ø8/8 Material: PEEK 	Pin To Rod Integrated Coupling REF: ZSSP04 Size: Ø5/6/8 Material: Aluminum 	5-hole Pin Clamp REF: 7100-0203 Size: Ø8/5-6 Material: Aluminum 
Tube To Rod Coupling REF: 7100-2001 Size: Ø20/8 Material: Aluminum 	Semi-Circular Curved Rod REF: 7100-0233/0213/0214 Size: Ø8x160/180/200 Material: Aluminum 	Straight Post REF: 7100-0205 Size: Ø8/0° Material: Stainless Steel 


COMBINATIONAL EXTERNAL FIXATOR SYSTEM

Parts & Instruments Of Large Fragment External Fixator Instrument Set

30° Post REF: 7100-0204 Size: Ø8/30° Material: Stainless Steel 	Connecting Rod REF: 7100-0234/7100-0210 Material: Aluminum/Carbon Fiber 
Compression/Distract Tube REF: 7100-2002 Size: Ø20 Material: Aluminum 	Compression/Distract Single Rod REF: 7100-2028 Size: Ø28 Material: Aluminum 
Compression/Distract Double Rods REF: 7100-2029 Size: Ø8 Material: Aluminum 	
Ankle Joint Movement Device REF: 7100-1903 Size: Ø8 Material: Stainless Steel 	Knee Joint Movement Device REF: 7100-1803 Size: Ø8 L/R Material: Stainless Steel 

Parts & Instruments Of Large Fragment External Fixator Instrument Set

Torque Wrench
REF: 7100-0216
Size: #7



Stabilization/Reduction Wrench
REF: 7100-0215
Size: #22




Wrench
REF: 7100-0230
Size: #5/6/7
Material: Aluminum




Hand Drill
REF: 7100-0219




T Wrench
REF: 7100-0220
Size: #7



T-Shaped Pin Wrench
REF: 7100-0221
Size: Ø5, Ø6



Drill Sleeve
REF: 7100-0209/0231
Size: Ø5/Ø6



Apex Pin Chuck
REF: 7100-0218
Size: Ø5, Ø6





Thumb Wheel
REF: 7100-0217
Size: Ø8/Ø7



Proximal Tibia External Fixator
REF: 7100-12
Used for: Proximal Tibia

REF	PRODUCT NAME	SPECIFICATION	QTY
7100-1201	Rod To Rod Coupling	Ø8/8	4
7100-1202	Pin To Rod Coupling	Ø8/5-6	3
7100-1203	5-hole Pin Clamp	Ø8/5-6	1
7100-1204	Post	Ø8/30°	2
7100-1205	Semi-Circular Curved Rod	Ø8 180	1
7100-1206	Connecting Rod	Ø8x250mm	2
7100-1207	Apex Pin Self-drilling	Ø5x150mm	5



Femoral Shaft External Fixator
REF: 7100-13
Used for: Femoral Shaft

REF	PRODUCT NAME	SPECIFICATION	QTY
7100-1301	Rod To Rod Coupling	Ø8/8	2
7100-1302	Pin To Rod Coupling	Ø8/5-6	4
7100-1303	Connecting Rod	Ø8x200mm	1
7100-1304	Connecting Rod	Ø8x250mm	2
7100-1305	Apex Pin Self-drilling	Ø5x150mm	4




Half Ring External Fixator
REF: 7100-14
Used for: Tibia, Perone

REF	PRODUCT NAME	SPECIFICATION	QTY
7100-1401	Rod To Rod Coupling	Ø8/8	6
7100-1402	5-hole Pin Clamp	Ø8/5-6	2
7100-1403	Post	Ø8/30°	4
7100-1404	Connecting Rod	Ø8x300mm	2
7100-1405	Connecting Rod	Ø8x250mm	1
7100-1406	Apex Pin Self-drilling	Ø5x150mm	4

Proximal Femur External Fixator

REF: 7100-15
Used for: Pertochanteric

REF	PRODUCT NAME	SPECIFICATION	QTY
7100-1501	Rod To Rod Coupling	Ø8/8	4
7100-1502	5-hole Pin Clamp	Ø8/5-6	2
7100-1503	Post	Ø8/30°	4
7100-1504	Connecting Rod	Ø8x300mm	2
7100-1505	Apex Pin Self-drilling	Ø5x150mm	2
7100-1506	Apex Pin Self-drilling	Ø5x180mm	2



Pelvic External Fixator


REF: 7100-16
Used for: Pelvic

REF	PRODUCT NAME	SPECIFICATION	QTY
7100-1601	Rod To Rod Coupling	Ø8/8	6
7100-1602	5-hole Pin Clamp	Ø8/5-6	2
7100-1603	Post	Ø8/30°	2
7100-1604	Connecting Rod	Ø8x300mm	2
7100-1605	Connecting Rod	Ø5x200mm	2
7100-1606	Apex Pin Self-drilling	Ø5x150mm	4

Ankle Joint External Fixator

REF: 7100-17
Used for: Ankle Joint

REF	PRODUCT NAME	SPECIFICATION	QTY
7100-1701	Rod To Rod Coupling	Ø8/8	3
7100-1702	Pin To Rod Coupling	Ø8/5-6	4
7100-1703	Connecting Rod	Ø8x250mm	2
7100-1704	Connecting Rod	Ø8x200mm	1
7100-1705	Apex Pin Self-drilling	Ø5x150mm	4



Knee Joint External Fixator

REF: 7100-18
Used for: Knee Joint

REF	PRODUCT NAME	SPECIFICATION	QTY
7100-1801	Rod To Rod Coupling	Ø8/8	2
7100-1802	Pin To Rod Coupling	Ø8/5-6	4
7100-1803	Knee Joint Movement Device	Ø8	1
7100-1804	5-hole Pin Clamp	Ø8/5-6	1
7100-1805	Post	Ø8/30°	1
7100-1806	Connecting Rod	Ø8x150mm	1
7100-1807	Apex Pin Self-drilling	Ø5x150mm	6



Ankle Fusion Joint External Fixator

REF: 7100-19
Used for: Ankle Joint


REF	PRODUCT NAME	SPECIFICATION	QTY
7100-1901	Pin To Rod Coupling	Ø5/3-4	1
7100-1902	Pin To Rod Coupling	Ø8/5-6	4
7100-1903	Ankle Joint Movement Device	Ø8	1
7100-1904	Apex Pin Self-drilling	Ø5x150mm	4
7100-1905	Apex Pin Self-drilling	Ø3x120mm	1



Micro-motion Type External Fixator

REF: 7100-20


REF	PRODUCT NAME	SPECIFICATION	QTY
7100-2001	Tube To Rod Coupling	Ø20/8	2
7100-2002	Compression/Distracton Tube	Ø20	1
7100-2003	5-hole Pin Clamp	Ø8/5-6	2
7100-2004	Post	Ø8/30°	2
7100-2005	Apex Pin Self-drilling	Ø5x150mm	4



Micro-motion Type External Fixator



REF: 7100-21

REF	PRODUCT NAME	SPECIFICATION	QTY
7100-2101	Tube To Rod Coupling	Ø20/8	2
7100-2102	Compression/Distracton Tube	Ø20	1
7100-2103	5-hole Pin Clamp	Ø8/5-6	2
7100-2104	Post	Ø8/30°	2
7100-2105	Apex Pin Self-drilling	Ø6x150mm	4





COMBINATIONAL EXTERNAL FIXATOR SYSTEM

ZKK External Fixator Instrument Set-I 8100-01

	REF	NO.	PRODUCT NAME	SPECIFICATION	QTY
	8100-0101	1	Ringent Rod To Rod Coupling	Ø8/8	4
	8100-0102	2	Ringent Rod To Rod Coupling	Ø8/4-5	8
	8100-0103	3	Ringent Rod To Rod Coupling	Ø11/11	4
	8100-0104	4	Ringent Rod To Rod Coupling	Ø11/5-6	12
	8100-0105	5	T Clamp	Ø11/5	2
	8100-0106	6	Shanz Pin	HA 4.0x90	4
	8100-0107			HA 5.0x110	4
	8100-0108			HA 5.0x130	6
	8100-0109			HA 6.0x130	6
	8100-0110			HA 5.0x150	2
	8100-0111			HA 6.0x150	2
	8100-0112	7	Drill Bits	Ø2.8x200	2
	8100-0113			Ø3.8x200	2
	8100-0114			Ø4.8x200	2
	8100-0115	8	Screw Driver	Ø4	1
	8100-0116			Ø5	1
	8100-0117	9	Tapered Trocar & Guide	Ø4	1
	8100-0118			Ø5	1
	8100-0119			Ø6	1
	8100-0120	10	Open Wrench	#10	1
	8100-0121	11	Pressurizer		1
	8100-0122	12	Connecting Rod	Ø8x150	2
	8100-0123			Ø8x200	2
	8100-0124			Ø8x250	2
	8100-0125			Ø8x300	2
	8100-0126			Ø11x150	2
	8100-0127			Ø11x200	2
	8100-0128			Ø11x250	2
	8100-0129			Ø11x300	2
	8100-0130			Ø11x350	2
	8100-0131	13	Aluminum Box		1


COMBINATIONAL EXTERNAL FIXATOR SYSTEM

ZAS External Fixator Instrument Set-II 8100-02

	REF	NO.	PRODUCT NAME	SPECIFICATION	QTY
	8100-0201	1	Pin To Rod Coupling	Ø11/5-6	12
	8100-0202	2	Ringent Rod To Rod Coupling	Ø11/11	9
	8100-0203	3	T Clamp	Ø11/5	2
	8100-0204	4	Shanz Pin	HA 5.0x110	4
	8100-0205			HA 5.0x130	6
	8100-0206			HA 5.0x150	2
	8100-0207			HA 6.0x130	6
	8100-0208			HA 6.0x150	2
	8100-0209	5	Drill Bits	Ø3.8x200	2
	8100-0210			Ø4.8x200	2
	8100-0211	6	Screw Driver	Ø5	1
	8100-0212			Ø6	1
	8100-0213	7	Tapered Trocar & Guide	Ø4	1
	8100-0214			Ø5	1
	8100-0215			Ø6	1
	8100-0216	8	Open Wrench	#10	2
	8100-0217	9	Pressurizer		1
	8100-0218	10	Connecting Rod	Ø11x150	4
	8100-0219			Ø11x200	4
	8100-0220			Ø11x250	4
	8100-0221			Ø11x300	4
	8100-0222			Ø11x350	4
	8100-0223	11	Aluminum Box		1

COMBINATIONAL EXTERNAL FIXATOR SYSTEM

ZKK External Fixator Instrument Set-III 8100-03



REF	NO.	PRODUCT NAME	SPECIFICATION	QTY
8100-0301	1	Pin To Rod Coupling	Ø8/3.5-4	8
8100-0302	2	Rod To Rod Coupling	Ø8/8	4
8100-0303	3	Rod To Rod Coupling	Ø11/11	4
8100-0304	4	Pin To Rod Coupling	Ø11/5-6	12
8100-0305	5	T Clamp	Ø11/5	2
8100-0306	6	Shanz Pin	HA 4.0x90	4
8100-0307			HA 5.0x110	4
8100-0308			HA 5.0x130	6
8100-0309			HA 6.0x130	6
8100-0310			HA 5.0x150	2
8100-0311	7	Drill Bits	HA 6.0x150	2
8100-0312			Ø2.8x200	2
8100-0313			Ø3.8x200	2
8100-0314	8	Screw Driver	Ø4.8x200	2
8100-0315			Ø4	1
8100-0316	9	Tapered Trocar & Guide	Ø5	1
8100-0317			Ø4	1
8100-0318			Ø5	1
8100-0319	10	Pressurizer	Ø6	1
8100-0320				
8100-0321				
8100-0322	11	Open Wrench	#10	2
8100-0323			Ø8x150	2
8100-0324			Ø8x200	2
8100-0325			Ø8x250	2
8100-0326			Ø8x300	2
8100-0327			Ø11x150	2
8100-0328			Ø11x200	2
8100-0329			Ø11x250	2
8100-0330			Ø11x300	2
8100-0331			Ø11x350	2
8100-0331	13	Aluminum Box		1

COMBINATIONAL EXTERNAL FIXATOR SYSTEM


Parts & Instrument Of ZKK, ZAS, ZIY External Fixator Instrument Set

Pin To Rod Coupling REF: 8100-0401 Size: Ø11/5 Material: Stainless Steel 	Pin To Rod Coupling REF: 8100-0402/0403 Size: Ø11/5-6, Ø8/3.5-4 Material: Aluminum 	T Clamp REF: 8100-0411 Size: Ø11/5 Material: Stainless Steel 
Rod To Rod Coupling REF: 8100-0414 Size: Ø11/11 Material: Stainless Steel 	Rod To Rod Coupling REF: 8100-0415/0416 Size: Ø11/11, Ø8/8 Material: Aluminum 	T Clamp REF: 8100-0412/0413 Size: Ø8/4, Ø11/5 Material: Aluminum 
Ringet Rod To Rod Coupling REF: 8100-0404/0405/0406 Size: Ø11/11, Ø11/8, Ø8/8 Material: Aluminum 	Ringet Pin To Rod Coupling REF: 8100-0407/0408 Size: Ø8/2.5-3.5, Ø8/4-5 Material: Aluminum 	Ringet Pin To Rod Coupling REF: 8100-0409/0410 Size: Ø11/5-6, Ø11/3-4 Material: Aluminum 

Parts & Instrument Of ZKK, ZAS, ZIY External Fixator Instrument Set

Integrated Pin To Rod Coupling REF: ZKTP01 Size: Ø8/3-4, Ø11/5-6 Material: Aluminum 	Integrated Rod To Rod Coupling REF: ZKTR01 Size: Ø8/8, Ø11/11 Material: Aluminum 	Ringent Rod To Rod Coupling REF: 8100-0419 Size: Ø11/11 Material: Aluminum 
Pin To Rod Coupling REF: 8100-0417/0418 Size: Ø8/3-4, Ø11/5-6 Material: Aluminum 	Semi-Circular Curved Rod REF: 8100-0420 Size: Ø8mm Material: Aluminum 	Trapezoidal Curved Rod REF: 8100-0422 Size: Ø8mm Material: Aluminum 
Connecting Rod REF: 8100-0423/0431 Size: 8x150/200/250/300 11x150/200/250/300/350 		

Parts & Instrument Of ZKK, ZAS, ZIY External Fixator Instrument Set

Curved Rod REF: 8100-0421 Size: Ø8mm Material: Aluminum 	Pressurizer REF: 8100-0432 
Compression/Distracton Rod REF: 8100-0433/0434 Size: Ø8 Material: Stainless Steel 	Open Spanner REF: ZZLGJ01 Size: 10# 
Tapered Trocar & Guide REF: ZZLGJ03 Size: Ø3mm, Ø4mm, Ø5mm, Ø6mm Note: Single or Triple is OK 	Hex Wrench REF: ZZLGJ04 Size: Ø2mm, Ø3mm, Ø4mm, Ø5mm, Ø6mm 
Screw Driver REF: ZZLGJ05 Size: Ø2mm, Ø2.5mm, Ø3mm, Ø3.5mm, Ø4mm, Ø5mm, Ø6mm 	

Half Ring Type External Fixator
REF: 8100-05
Used for: Tibia & perone open & closed fracture

REF	PRODUCT NAME	SPECIFICATION	QTY
8100-0501	Ringent Rod To Rod Coupling	Ø8/8	6
8100-0502	Ringent Pin To Rod Coupling	Ø8/4-5	6
8100-0503	Semi-Circular Curved Rod	Ø8	2
8100-0504	Connecting Rod	Ø8x300	1
8100-0505	Compression/Distraction Rod	Ø8	2
8100-0506	Shanz Pin	HA 4.0x120	6



Z Frame(T Type) External Fixator
REF: 8100-08
Used for: Femur comminuted fracture

REF	PRODUCT NAME	SPECIFICATION	QTY
8100-0801	Ringent Rod To Rod Coupling	Ø8/8	2
8100-0802	Ringent Pin To Rod Coupling	Ø8/4-5	4
8100-0803	T Clamp	Ø8/4	1
8100-0804	Connecting Rod	Ø8x100	1
8100-0805	Connecting Rod	Ø8x150	2
8100-0806	Shanz Pin	HA 4.0x120	6



Condylar External Fixator
REF: 8100-06
Used for: Humeral epicondyle, tibial epicondyle, distal fracture

REF	PRODUCT NAME	SPECIFICATION	QTY
8100-0601	Ringent Rod To Rod Coupling	Ø8/8	3
8100-0602	Ringent Pin To Rod Coupling	Ø8/4-5	6
8100-0603	Semi-Circular Curved Rod	Ø8	1
8100-0604	Trapezoidal Curved Rod	Ø8	1
8100-0605	Connecting Rod	Ø8x250	1
8100-0606	Shanz Pin	HA 4.0x120	6

T Frame External Fixator
REF: 8100-09
Used for: Femur, Tibia

REF	PRODUCT NAME	SPECIFICATION	QTY
8100-0901	Pin To Rod Coupling	Ø11/5-6	4
8100-0902	T Clamp	Ø11/5	1
8100-0903	Connecting Rod	Ø11x300	1
8100-0904	Shanz Pin	Ø5x130	4
8100-0905	Shanz Pin	HBØ5x150	2



Straight Type External Fixator
REF: 8100-10
Used for: Femur, Tibia

REF	PRODUCT NAME	SPECIFICATION	QTY
8100-1001	Ringent Rod To Rod Coupling	Ø11/11	1
8100-1002	Pin To Rod Coupling	Ø11/5-6	4
8100-1003	Connecting Rod	Ø11x200	2
8100-1004	Shanz Pin	Ø5x130	4




Tibial Joint Side External Fixator
REF: 8100-07
Used for: Distal tibial joint, proximal fracture

REF	PRODUCT NAME	SPECIFICATION	QTY
8100-0701	Ringent Rod To Rod Coupling	Ø8/8	3
8100-0702	Ringent Pin To Rod Coupling	Ø8/4-5	7
8100-0703	Semi-Circular Curved Rod	Ø8	1
8100-0704	Curved Rod	Ø8	1
8100-0705	Connecting Rod	Ø8x300	1
8100-0706	Shanz Pin	HA 4.0x120	7



Z Type External Fixator
REF: 8100-11
Used for: Femur, Tibia

REF	PRODUCT NAME	SPECIFICATION	QTY
8100-1101	Rod To Rod Coupling	Ø11/11	2
8100-1102	Pin To Rod Coupling	Ø11/5-6	4
8100-1103	Connecting Rod	Ø11x200	2
8100-1104	Connecting Rod	Ø11x150	1
8100-1105	Shanz Pin	Ø5x130	4



» COMBINATIONAL EXTERNAL FIXATOR SYSTEM


T Type External Fixator
REF: 8100-12
Used for: Tibia comminuted fracture

REF	PRODUCT NAME	SPECIFICATION	QTY
8100-1201	Pin To Rod Coupling	Ø11/5-6	4
8100-1202	T Clamp	Ø11/5	1
8100-1203	Connecting Rod	Ø11x300	1
8100-1204	Shanz Pin	Ø5x130	4
8100-1205	Shanz Pin	HAØ5x150	2



Z Type(Enhancement Mode) External Fixator
REF: 8100-13
Used for: Femur, Tibia

REF	PRODUCT NAME	SPECIFICATION	QTY
8100-1301	Rod To Rod Coupling	Ø11/11	2
8100-1302	Pin To Rod Coupling	Ø11/5-6	4
8100-1303	Connecting Rod	Ø11x200	2
8100-1304	Connecting Rod	Ø11x150	1
8100-1305	Shanz Pin	Ø5x130	4



Straight Type External Fixator
REF: 8100-14
Used for: Tibia

REF	PRODUCT NAME	SPECIFICATION	QTY
8100-1401	Pin To Rod Coupling	Ø11/5-6	4
8100-1402	Connecting Rod	Ø11x250	1
8100-1403	Shanz Pin	Ø5x130	4



Pelvic Fixator External Fixator
REF: 8100-15
Used for: Pelvic


REF	PRODUCT NAME	SPECIFICATION	QTY
8100-1501	Ringent Rod To Rod Coupling	Ø11/11	2
8100-1502	Ringent Pin To Rod Coupling	Ø11/5-6	4
8100-1503	Ringent Rod To Rod Coupling	Ø11/8	2
8100-1504	Connecting Rod	Ø11x200	2
8100-1505	Connecting Rod	Ø11x300	1
8100-1506	Compression/Disraction Rod	Ø8	1
8100-1507	Shanz Pin	Ø5x150	4

» COMBINATIONAL EXTERNAL FIXATOR SYSTEM

Universal Pin/Rod To Rod Coupling

Z Type External Fixator
REF: ZWX001

REF	PRODUCT NAME	SPECIFICATION	QTY
ZWXR01	Universal Rod To Rod Coupling	Ø11/11	2
ZWXP01	Universal Pin To Rod Coupling	Ø11/6	4
ZIYL01	Connecting Rod	Ø11x200	2
ZIYL01	Connecting Rod	Ø11x150	1
S\$1060134S	Shanz Pin	Ø6x130	4



Universal Rod To Rod Coupling
REF: ZWXR01
Size: Ø11/11, Ø8/8




Universal Pin To Rod Coupling
REF: ZWXP01
Size: Ø11/6, Ø8/4



Carbon Fiber Composite Ring

Three-Dimensional Type External Fixator
REF: ZHH001

REF	PRODUCT NAME	SPECIFICATION	QTY
ZHHL01	Rings	L	1
ZIYL01	Connecting Rod	Ø11/250	2
ZHHR01	Rod To Ring Coupling	Ø11/11	2
ZHHP01	Pin To Ring Coupling	Ø2.0	4
ZKKP01	Ringent Pin To Rod Coupling	Ø11/5-6	2
ZKKR01	Ringent Rod To Rod Coupling	Ø11/11	1
SS10603345	Threaded Rod	Ø6x130	2
SSS020250	K-Wire	Ø2.0x250	2




Tight Wire Clamp
REF: ZHH-G01



Rings
REF: ZHHL01



Rod To Ring Coupling
REF: ZHHR01
Size: Ø11/11



Pin To Ring Coupling
REF: ZHHP01
Size: Ø2.4



- Straight Type External Fixator
- L
REF: 9100-0101
Matching Bone Screw: Ø6.0x130mm 4 pcs
Matching Instruments: 6mm hex wrench, 6mm screw driver
- M
REF: 9100-0102
Matching Bone Screw: Ø6.0x120mm 4 pcs
Matching Instruments: 6mm hex wrench, 6mm screw driver
(Select 5mm bone screw)
- S
REF: 9100-0103
Matching Bone Screw: Ø5.0x110mm 4 pcs
Matching Instruments: 6mm hex wrench, 5mm screw driver
- XS
REF: 9100-0104
Matching Bone Screw: Ø4.0x90mm 4 pcs
Matching Instruments: 4mm hex wrench, 4mm screw driver



- T Type External Fixator
- L
REF: 9100-0201
Matching Bone Screw: Ø6.0x130mm 2 pcs, HBØ6.0x150mm 2 pcs
Matching Instruments: 6mm hex wrench, 6mm screw driver
- M
REF: 9100-0202
Matching Bone Screw: Ø6.0x120mm 4 pcs, HBØ6.0x130mm 2 pcs
Matching Instruments: 6mm hex wrench, 6mm screw driver
- S
REF: 9100-0203
Matching Bone Screw: Ø5.0x110mm 2 pcs, HBØ5.0x120mm 2 pcs
Matching Instruments: 5mm hex wrench, 5mm screw driver
- XS
REF: 9100-0204
Matching Bone Screw: Ø4.0x90mm 4 pcs
Matching Instruments: 4mm hex wrench, 4mm screw driver



- Ankle Joint External Fixator
- L
REF: 9100-0301
Matching Bone Screw: Ø6.0x130mm 4 pcs
Matching Instruments: 6mm hex wrench, 6mm screw driver
- M
REF: 9100-0302
Matching Bone Screw: Ø5.0x130mm 4 pcs
Matching Instruments: 6mm hex wrench, 5mm screw driver



Straight Type With Middle Pins External Fixator

L
REF: 9100-0401
Matching Bone Screw: Ø6.0x130mm 6 pcs
Matching Instruments: 6mm hex wrench, 6mm screw driver

M
REF: 9100-0402
Matching Bone Screw: Ø6.0x120mm 6 pcs
Matching Instruments: 6mm hex wrench, 6mm screw driver
(Select 5mm bone screw)

S
REF: 9100-0403
Matching Bone Screw: Ø5.0x110mm 6 pcs
Matching Instruments: 6mm hex wrench, 5mm screw driver



T Type With Middle Pins External Fixator

L
REF: 9100-0501
Matching Bone Screw: Ø6.0x130mm 4 pcs, HBØ6.0x150mm 2 pcs
Matching Instruments: 6mm hex wrench, 6mm screw driver

M
REF: 9100-0502
Matching Bone Screw: Ø6.0x120mm 4 pcs, HBØ6.0x130mm 2 pcs
Matching Instruments: 6mm hex wrench, 6mm screw driver

S
REF: 9100-0503
Matching Bone Screw: Ø5.0x110mm 4 pcs, HBØ5.0x120mm 2 pcs
Matching Instruments: 5mm hex wrench, 5mm screw driver



3D Type External Fixator

L
REF: 9100-0601
Matching Bone Screw: Ø6.0x130mm 4 pcs
Matching Instruments: 10-14mm open spanner, 6mm hex wrench

M
REF: 9100-0602
Matching Bone Screw: Ø5.0x120mm 4 pcs
Matching Instruments: 10-14mm open spanner, 5mm hex wrench

S
REF: 9100-0603
Matching Bone Screw: Ø5.0x110mm 4 pcs
Matching Instruments: 10-14mm open spanner, 5mm hex wrench

XS
REF: 9100-0604
Matching Bone Screw: Ø4.0x90mm 4 pcs
Matching Instruments: 10mm open spanner, 4mm hex wrench



Three-Selection Straight Type External Fixator

REF: 9100-0701
Matching Bone Screw: Ø6.0x130mm 6 pcs
Matching Instruments: 6mm hex wrench, 6mm screw driver



Three-Selection T-Shaped Type External Fixator

REF: 9100-0801
Matching Bone Screw: Ø6.0x130mm 4 pcs, HBØ6.0x150mm 2 pcs
Matching Instruments: 6mm hex wrench, 6mm screw driver



Joint Type External Fixator

REF: 9100-0901
Matching Bone Screw: Ø6.0x130mm 2 pcs, HBØ6.0x180mm 2 pcs
Matching Instruments: 6mm hex wrench, 6mm screw driver




Oblique Type External Fixator

REF: 9100-1001
Matching Bone Screw: Ø6.0x150mm 2 pcs, HBØ6.0x180mm 2 pcs
Matching Instruments: 6mm hex wrench, 6mm screw driver



» COMMON EXTERNAL FIXATOR SYSTEM


Dynamic External Fixator
REF: 9100-1101
Matching Bone Screw: Ø4.0x250mm 2 pcs, HBØ5.0x250mm 3 pcs
Matching Instruments: 10mm open spanner, 4&5 mm screw driver



Stainless Steel Universal Type External Fixator-I
L
REF: 9100-1201
Matching Bone Screw: Ø6.0x130mm 4 pcs
Matching Instruments: 13-16mm open spanner, 4&5 mm hex wrench 6mm screw driver

M
REF: 9100-1202
Matching Bone Screw: Ø5.0x120mm 4 pcs
Matching Instruments: 10-14mm open spanner, 4&5 mm hex wrench 5mm screw driver


S
REF: 9100-1203
Matching Bone Screw: Ø5.0x110mm 4 pcs
Matching Instruments: 10-14mm open spanner, 4&5 mm hex wrench 5mm screw driver



Stainless Steel Straight Type External Fixator-II
L
REF: 9100-1301
Matching Bone Screw: Ø6.0x130mm 4 pcs
Matching Instruments: 10-16mm open spanner, 4mm hex wrench 6mm screw driver


M
REF: 9100-1302
Matching Bone Screw: Ø5.0x120mm 4 pcs
Matching Instruments: 10-14mm open spanner, 4mm hex wrench 5mm screw driver

S
REF: 9100-1303
Matching Bone Screw: Ø5.0x110mm 4 pcs
Matching Instruments: 10-14mm open spanner, 4mm hex wrench 5mm screw driver




» COMMON EXTERNAL FIXATOR SYSTEM


Mini-External Fixator (Finger & Palm)
REF: 9100-1401
Used for: Finger, Palm



Mini-External Fixator (Finger & Palm)
REF: 9100-1402
Used for: Finger, Palm




Wrist Joint External Fixator
REF: 9100-1501
Used for: Wrist Joint




Semi-Circular Type (Standard) External Fixator
REF: 9100-1601
Used for: Tibia&perone open&closed fracture

REF	PRODUCT NAME	SPECIFICATION	QTY
ZZLR01	Rod To Rod Coupling	Ø8/8	6
ZZLP01	Pin To Rod Coupling	Ø8/4-5	6
ZZLL02	Semi-Circular Curved Rod	Ø8mm L	2
ZZLL01	Rod	Ø8x300	1
ZZLL04	Compression/Distracton Rod	Ø8	2
SS10401240	Shanz Pin	Ø4x120	2
SS5040250	K-wire	Ø4x250	2



Pelvic Fixator
REF: 9100-1701
Used for: Pelvic


REF	PRODUCT NAME	SPECIFICATION	QTY
ZZLR01	Rod To Rod Coupling	Ø8/8	6
ZZLP01	Pin To Rod Coupling	Ø8/4-5	6
ZZLL03	Semi-Circular Curved Rod	Ø8mm	2
ZZLL04	Compression/Distracton Rod	Ø8	1
ZZLL01	Rod	Ø8x180	2
SS10501345	Shanz Pin	Ø5x130	6



Double Rods Tuberosity External Fixator

REF: 9100-1801
Used for: Pertrochanteric

REF	PRODUCT NAME	SPECIFICATION	QTY
ZZLP02	Double Rods Pin Clamp	Ø8/4-5	5
ZZLL01	Rod	Ø8x300	2
SS10501345	Shanz Pin	Ø5x130	2
SS10501850	Shanz Pin	Ø5x180	3



Double Rod Single Side External Fixator

REF: 9100-1901
Used for: Femur, Tibia

REF	PRODUCT NAME	SPECIFICATION	QTY
ZZLP02	Pin To Double Rod Coupling	Ø8/4-5	4
ZZLL01	Rod	Ø8x250	2
SS10501345	Shanz Pin	Ø5x130	4



Humeral Lateral Condylar Neck External Fixator

REF: 9100-2001
Used for: Humeral lateral condylar neck, proximal & distal humer


REF	PRODUCT NAME	SPECIFICATION	QTY
ZZLP01	Pin To Rod Coupling	Ø8/4-5	4
ZZLR01	Rod To Rod Coupling	Ø8/8	1
ZZLL01	Rod	Ø8x200	1
ZZLL01	Rod	Ø8x150	1
SS10501345	Shanz Pin	Ø5x130	4



Humeral Frame External Fixator

REF: 9100-2101
Used for: Humer


REF	PRODUCT NAME	SPECIFICATION	QTY
ZZLR01	Rod To Rod Coupling	Ø8/8	2
ZZLP01	Pin To Rod Coupling	Ø8/4-5	6
ZZLL01	Rod	Ø8x250	2
ZZLL01	Rod	Ø8x200	1
SS10401240	Shanz Pin	Ø4x120	6



Femoral Condylar External Fixator

REF: 9100-2201
Used for: Femoral supracondylar& intercondylar fracture


REF	PRODUCT NAME	SPECIFICATION	QTY
ZZLR01	Rod To Rod Coupling	Ø8/8	3
ZZLP02	Rod To Double Rod Coupling	Ø8/4-5	2
ZZLP01	Pin To Rod Coupling	Ø8/4-5	5
ZZLL02	Semi-Circular Curved Rod	Ø8mm	1
ZZLL01	Rod	Ø8x300	2
ZZLL03	Curved Rod	Ø8mm Long	1
SS10501345	Shanz Pin	Ø5x130	7



Femoral Shaft External Fixator

REF: 9100-2301
Used for: Femoral shaft 1/3 to the joint comminuted fracture

REF	PRODUCT NAME	SPECIFICATION	QTY
ZZLR01	Rod To Rod Coupling	Ø8/8	3
ZZLP02	Pin To Double Rod Coupling	Ø8/4-5	4
ZZLP01	Pin To Rod Coupling	Ø8/4-5	3
ZZLL03	Semi-Circular Curved Rod	Ø8mm Long	2
ZZLL01	Rod	Ø8x300	2
SS10501345	Shanz Pin	Ø5x130	4
SS 5040250	K-wire	Ø4x250	1



Ankle Joint External Fixator

REF: 9100-2401
Used for: Distal tibia, Ankle joint



REF	PRODUCT NAME	SPECIFICATION	QTY
ZZLR01	Rod To Rod Coupling	Ø8/8	6
ZZLP01	Pin To Rod Coupling	Ø8-4/5	8
ZZLL04	Compression/Distracton Rod	Ø8	2
ZZLL01	Rod	Ø8x350	2
ZZLL01	Rod	Ø8x250	2
SS10401240	Shanz Pin	Ø4x120	4
SS 5040250	K-wire	Ø4x250	2



Tibial Plateau Frame External Fixator

REF: 9100-2501
Used for: Proximal tibia


REF	PRODUCT NAME	SPECIFICATION	QTY
ZZLR01	Rod To Rod Coupling	Ø8/8	3
ZZLP01	Pin To Rod Coupling	Ø8/4-5	8
ZZLL01	Rod	Ø8x300	1
ZZLL02	Semi-Circular Curved Rod	Ø8mm L	1
ZZLL03	Curved Rod	Ø8mm Long	1
SS10501345	Shanz Pin	Ø5x130	6
SS 5040250	K-wire	Ø4x250	1



Humeral Condylar External Fixator

REF: 9100-2601
Used for: Humeral&tibial lateral epicondyle

REF	PRODUCT NAME	SPECIFICATION	QTY
ZZLR01	Rod To Rod Coupling	Ø8/8	3
ZZLP01	Pin To Rod Coupling	Ø8/4-5	5
ZZLL01	Rod	Ø8x250	1
ZZLL02	Semi-Circular Curved Rod	Ø8mm L	1
ZZLL03	Curved Rod	Ø8mm Short	1
SS10401240	Shanz Pin	Ø4x120	3
SS 5040250	K-wire	Ø4x250	1



Knee Joint External Fixator

REF: 9100-2701
Used for: Spanning joint,
distal femur&proximal tibia

REF	PRODUCT NAME	SPECIFICATION	QTY
ZZLR01	Rod To Rod Coupling	Ø8/8	8
ZZLP01	Pin To Rod Coupling	Ø8/4-5	14
ZZLL02	Semi-Circular Curved Rod	Ø8mm L	2
ZZLL03	Curved Rod	Ø8mm Long	2
ZZLL06	Spaning Joint Connecting Rod	Ø8	2
SS10501345	Shanz Pin	Ø5x130	14

Parts And Tools For Combination External Fixator ZZL

Single Rod Clamp

REF: ZZLP01
Size: 6Ø/2.5-3.5mm, Ø8/3-4mm,
Ø8/4-5mm, Ø8/5-6mm



Double Rods Clamp

REF: ZZLR01
Size: 6Ø/6mm, Ø8/8mm



Double Rods Pin Clamp

REF: ZZLP02
Size: 8Ø/3-4mm, Ø8/4-5mm, Ø8/5-6mm




Double Rod To Rod Clamp

REF: ZZLR02
Size: 6Ø/6mm, Ø8/8mm



Semi-Circular Curved Connecting Rod

REF: ZZLL02
Size: Ø6mm L, Ø6mm S,
Ø8mm L, Ø8mm S
Material: Stainless Steel



Trapezoidal Curved Connecting Rod

REF: ZZLL05
Size: Ø6mm, Ø8mm
Material: Stainless Steel



Curved Connecting Rod

REF: ZZLL03
Size: Ø6mm L, Ø6mm S, Ø8mm L, Ø8mm S
Material: Stainless Steel



Rod

REF: ZZLL01



» COMMON EXTERNAL FIXATOR SYSTEM

Parts And Tools For Combination External Fixator ZZL

Compression/Distracton Rod
REF: ZZLL04
Size: Ø8
Material: Stainless Steel



Spanning Joint Connecting Rod
REF: ZZLL06
Size: Ø8
Material: Stainless Steel



Open Spanner
REF: ZZLGJ01
Size: 10#



Open Spanner
REF: ZZLGJ02
Size: 10-13




Tapered Trocar & Guide
REF: ZZLGJ03
Size: Ø3mm, Ø4mm, Ø5mm, Ø6mm
Note: Single or Triple is OK



Hex Wrench
REF: ZZLGJ04
Size: Ø2mm, Ø3mm, Ø4mm, Ø5mm, Ø6mm



Screw Driver
REF: ZZLGJ05
Size: Ø2mm, Ø2.5mm, Ø3mm, Ø3.5mm, Ø4mm, Ø5mm, Ø6mm




» ACCESSORIES OF EXTERNAL FIXATOR SYSTEM


Shanz Pin
REF: 1001-0101/0184
Diameter: 2.0/2.5/3.0/3.5/4.0/5.0/6.0mm
Length: 50-250mm



Apex Pin Self-Drilling
REF: 1001-0201/0325
Diameter: 3.0/4.0/5.0/6.0mm
Length: 60-250mm



Threaded Guide Wire
REF: 1001-0301/0325
Diameter: 1.2/1.5/2.0/2.5/3.0mm
Length: 100/120/150/250/300mm



Wire With Olive
REF: 1001-0401/0402/0403
Specification: 1.5x350/2.0x350/2.5x350mm




» ACCESSORIES OF EXTERNAL FIXATOR SYSTEM

K-Wire
REF: 1001-0501/0536
Diameter: 0.8/1.0/1.2/1.5/1.6/1.8/2.0/2.5/3.0/3.5/4.0/4.5mm
Length: 100/120/150/250/300mm



M6 Threaded Rod
REF: 1001-0601/0613
Length: 40-400mm



Drill Bit
REF: 1001-0701/0734
Diameter: 1.0/1.2/1.5/1.8/2.0/2.5/2.8/3.0/3.2/3.5/4.0/4.2/4.8/5.0mm
Length: 115-200mm



Connecting Rod
REF: 1001-0801/0827
Diameter: 3/5/6/8/11mm
Length: 40-400mm
Material: Carbon fiber/Aluminum/Stainless Steel





Orthopromed Inc.

Address: 964 E. Badillo Street #2140
Covina CA 91724 **UNITED STATES**

Phone: 302-613-4333

Web: www.orthopromed.com

Email: info@orthopromed.com

After Sales: support@orthopromed.com