




Orthopromed Inc.



Veterinary Orthopedics

www.orthopromed.com



CONTENTS

Products	Page
SOP / TPL System	3
Locking Plate System	8
Non Locking Plate System	18
Orthopedic Instruments	34
External Fixation System	43
Power Tool System	46
General Surgery Instruments	49
Dental Instruments	53

» SOP / TPL SYSTEM

2.0mm SOP System



Product Code	Specification	Material
GA0012020	20 Holes	Stainless Steel
Note: Match with 2.0 locking screw		

2.4mm SOP System



Product Code	Specification	Material
GA0012420	20 Holes	Stainless Steel
Note: Match with 2.4 locking screw		

2.7mm SOP System



Product Code	Specification	Material
GA0012720	20 Holes	Stainless Steel
Note: Match with 2.7 locking screw		

TPLO 2.0 Curved Plate



Product Code	Length / Thickness(mm)	Specification	Material
GA0022001	30/2.3	3+3 Holes L	Pure Titanium
GA0022002	30/2.3	3+3 Holes R	
Note: Match with 2.0 locking screw and 2.0 cortical screw, photo for the left			

TPLO 2.4 Curved Plate



Product Code	Length / Thickness(mm)	Specification	Material
GA0022401	35/2.3	3+3 Holes L	Pure Titanium
GA0022402	35/2.3	3+3 Holes R	
Note: Match with 2.4 locking screw and 2.4 cortical screw, photo for the right			

TPLO 2.7 Curved Plate



Product Code	Length / Thickness(mm)	Specification	Material
GA0022701	46/2.5	3+3 Holes L	Pure Titanium
GA0022702	46/2.5	3+3 Holes R	
Note: Match with 2.7 locking screw and 2.7 cortical screw, photo for the left			

» SOP / TPL SYSTEM

TPLO 3.5 Curved Plate



Product Code	Length / Thickness(mm)	Specification	Material
GA0023501	64/3.7	3+3 Holes L	Pure Titanium
GA0023502	64/3.7	3+3 Holes R	

Note: Match with 3.5 locking screw and 3.5 cortical screw, photo for the right

Arthodesis Of Wrist Joint Locking Plate (Straight Type)



Product Code	Length(mm)	Material
GA0032001	68	Pure Titanium

Note: Match with 2.0/2.4 locking screw and 2.0/2.4 cortical screw

Arthodesis Of Wrist Joint Locking Plate (Straight Type)



Product Code	Length(mm)	Material
GA0032002	82	Pure Titanium

Note: Match with 2.0/2.7 locking screw and 2.0/2.7 cortical screw

Arthodesis Of Tarsal Joint Locking Plate (L Type)



Product Code	Length(mm)	Material
GA0042001	135° L	Pure Titanium
GA0042002	135° R	

Note: Match with 2.0 locking screw and 2.0 cortical screw, photo for the left

Arthodesis Of Tarsal Joint Locking Plate (L Type)



Product Code	Length(mm)	Material
GA0042701	135° L	Pure Titanium
GA0042702	135° R	

Note: Match with 2.4/2.7 locking screw and 2.4/2.7 cortical screw, photo for the right

TType Locking Plate 2.0



Product Code	Length/Width(mm)	Material
GA00520	50/5.3	Pure Titanium

Note: Match with 2.0 locking screw and 2.0 cortical screw

TType Locking Plate 2.4/2.7



Product Code	Length/Width(mm)	Material
GA0052427	55/7	Pure Titanium

Note: Match with 2.4/2.7 locking screw and 2.4/2.7 cortical screw

» SOP / TPL SYSTEM

2.0 Locking Screw



Product Code	Specification(mm)	Material
GB0012006-GB0012026	2.0x8/10/12/14/16/18/20/22/24/26	Titanium Alloy

Note: Match T6 toxy screw driver

2.4 Locking Screw



Product Code	Specification(mm)	Material
GB0012408-GB0012430	2.4x8/10/12/14/16/18/20/22/24/26/28/30	Titanium Alloy

Note: Match T8 toxy screw driver

2.7 Locking Screw



Product Code	Specification(mm)	Material
GB0012710-GB0012734	2.7x10/12/14/16/18/20/22/24/26/28/30/32/34	Titanium Alloy

Note: Match T8 toxy screw driver

3.5 Locking Screw



Product Code	Specification(mm)	Material
GB0013510-GB0013560	3.5x10/12/14/16/18/20/22/24/26/28/30/32/34/36/38/40/45/50/55/60	Titanium Alloy

Note: Match T15 toxy screw driver

2.0 Cortical Self-Tapping Screw



Product Code	Specification(mm)	Material
GB0012006-GC0012026	2.0x6/7/8/10/12/14/16/18/20/22/24/26	Titanium Alloy

Note: Match T6 toxy screw driver

2.4 Cortical Self-Tapping Screw



Product Code	Specification(mm)	Material
GC0012410-GC0012432	2.4x10/12/14/16/18/20/22/24/26/28/30/32	Titanium Alloy
Note: Match T6 toxy screw driver		

2.7 Cortical Self-Tapping Screw



Product Code	Specification(mm)	Material
GC0012712-GC0012734	2.7x12/14/16/18/20/22/24/26/28/30/32/34	Titanium Alloy
Note: Match T8 toxy screw driver		

3.5 Cortical Self-Tapping Screw



Product Code	Specification(mm)	Material
GC0013510-GC0013560	3.5x10/12/14/16/18/20/22/24/26/28/30/32/34/36/38/40/45/50/55/60	Titanium Alloy
Note: Match T15 toxy screw driver		

Locking Drill Guide




Product Code	Specification(mm)
GD00120	2.0mm
GD00124	2.4mm
GD00127	2.7mm
GD00135	3.5mm

Quick Coupling Screw Driver




Product Code	Specification(mm)
GE001	T6
GE002	T8
GE003	T15

2.0/2.4 TPLO Insrument Kit



NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION	QTY
1	A008-01	Drill Bit (AO)	Ø1.5	1
2	A008-02	Drill Bit (AO)	Ø1.8	1
3	A008-03	Locking Drill Guide	1.5mm	1
4	A008-04	Locking Drill Guide	1.5mm	1
5	A008-05	Double Drill Guide	1.5/2.0	1
6	A008-06	Double Drill Guide	1.8/2.4	1
7	A008-07	Torx Head Screw Driver	T6	1
8	A008-08	Torx Head Screw Driver	T8	1
9	A008-09	AO Quick Coupling Handle		1
10	A008-10	Torque Limited AO Coupling Handle	0.4Nm	1
11	A008-11	Depth Gauge	30mm	1
12	A008-12	Mini TPLO Target Device	Small	1
13	A008-13	Proltractor	Small	1
14	A008-14	Screw Sterilization Box		1
15	A008-15	TPL Plate	2.0mm	2
16	A008-16	TPL Plate	2.4mm	2
17	A008-17	Tool Box		1
18	A008-18	Poritioning Pin	2.5x120mm	2


2.7/3.5 TPLO Insrument Kit



NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION	QTY
1	A009-01	Drill Bit (AO)	Ø2.0	1
2	A009-02	Drill Bit (AO)	Ø2.5	1
3	A009-03	Drill Bit (AO)	Ø2.7	1
4	A009-04	Locking Drill Guide	2.0mm	1
5	A009-05	Locking Drill Guide	2.7mm	1
6	A009-06	Double Drill Guide	2.0/2.7	1
7	A009-07	Double Drill Guide	2.5/3.5	1
8	A009-08	Torx Head Screw Driver	T8	1
9	A009-09	Torx Head Screw Driver	T15	1
10	A009-10	AO Quick Coupling Handle		1
11	A009-11	Torque Limited AO Coupling Handle	0.8NM	1
12	A009-12	Torque Handle	1.5Nm	
13	A009-13	Depth Gauge	60mm	1
14	A009-14	TPLO Target Device	Large	1
15	A009-15	Proltractor		1
16	A009-16	Screw Sterilization Box		1
17	A009-17	TPL Plate	2.7mm	2
18	A009-18	TPL Plate	3.5mm	2
19	A009-19	Poritioning Pin	2.5x120mm	2
20	A009-20	Tool Box		1

» LOCKING PLATE SYSTEM

1.6mm Reconstruction Locking Plate AA00116



Product Code	Thickness/WIDTH Length (mm)	Specification	Material
AA0011630	1.2/3.5/135	30 Holes	Pure Titanium
AA0011626	1.5/4.2/130	26 Holes	

Note: Used with 1.6 locking screw and 1.1 cortical screw


2.4mm Reconstruction Locking Plate AA00124



Product Code	Thickness/WIDTH Length (mm)	Specification	Material
AA0012443	1.8/5/236.5	43 Holes	Pure Titanium
AA0012434	2/6.5/238	34 Holes	

Note: Used with 2.4 locking screw and 1.5 cortical screw


3.2mm Reconstruction Locking Plate AA00132



Product Code	Thickness/WIDTH Length (mm)	Specification	Material
AA0013226	2.5/7.5/234	26 Holes	Pure Titanium
AA0013222	3/8.5/220	22 Holes	

Note: Used with 3.2 locking screw and 2.4 cortical screw

4.0mm Reconstruction Locking Plate AA00140



Product Code	Thickness/WIDTH Length (mm)	Specification	Material
AA0014022	3.5/10/240	20 Holes	Pure Titanium
AA0014018	4/11.5/234	18 Holes	

Note: Used with 4.0 locking screw and 2.7 cortical screw

» LOCKING PLATE SYSTEM


1.6mm T Type Locking Plate AA00216



Product Code	HEADWIDTH/WIDTH/Length (mm)	Specification	Material
AA0021601	8/3.5/30	2+5 Holes	Pure Titanium
AA0021602	9/4.2/33	2+5 Holes	
AA0021603	7/4.2/34	3+5 Holes	
AA0021604	6.5/3.5/31	3+5 Holes	

Note: Used with 1.6 locking screw and 1.1 cortical screw

2.4mm T Type Locking Plate AA00224



Product Code	HEADWIDTH/WIDTH/Length (mm)	Specification	Material
AA0022401	11/5/38	2+5 Holes	Pure Titanium
AA0022402	13/6.5/48	2+5 Holes	
AA0022403	9/5/39	3+5 Holes	

Note: Used with 2.4 locking screw and 1.5 cortical screw


3.2mm T Type Locking Plate AA00232



Product Code	HEADWIDTH/WIDTH/Length (mm)	Specification	Material
AA0023201	16/7.5/59	2+5 Holes	Pure Titanium
AA0023202	19/8.5/68	2+5 Holes	

Note: Used with 3.2 locking screw and 2.4 cortical screw

1.1mm T Type Plate(A Series) AA00311



Product Code	HEADWIDTH/WIDTH/Length (mm)	Specification	Material
AA0031101	1.0/3/16	2+3 Holes	Pure Titanium
AA0031102	1.0/3/20	2+4 Holes	
AA0031103	1.0/3/24	2+5 Holes	

Note: A series plate head to shaft distance is 4.6mm, used with 1.1mm cortical screw

» LOCKING PLATE SYSTEM

1.1mm T Type Plate(B Series) AA00411



Product Code	HEADWIDTH/WIDTH/Length (mm)	Specification	Material
AA0041101	1.0/3/18	2+3 Holes	Pure Titanium
AA0041102	1.0/3/22	2+4 Holes	
AA0041103	1.0/3/26	2+5 Holes	

Note: B series plate head to shaft distance is 6.6mm, Used with 1.1 cortical screw


1.1mm T Type Plate(C Series) AA00511



Product Code	HEADWIDTH/WIDTH/Length (mm)	Specification	Material
AA0051101	1.0/3/21	3+3 Holes	Pure Titanium
AA0051102	1.0/3/27	3+4 Holes	

Note: C series plate head to shaft distance is 7.4/9.4mm, Used with 1.1 cortical screw

1.1mm Cuttable Plate AA00611



Product Code	HEADWIDTH/WIDTH/Length (mm)	Specification	Material
AA0061101	1.0/3/90	30 Holes	Pure Titanium
AA0061102	1.0/3/91	23 Holes	

Note: Used with 1.1 cortical screw


1.1mm Reconstruction Cuttable Plate AA00711



Product Code	HEADWIDTH/WIDTH/Length (mm)	Specification	Material
AA0071101	1.0/3/90	30 Holes	Pure Titanium

Note: Used with 1.1 cortical screw

1.1mm Neutral Cuttable Plate AA00811



Product Code	HEADWIDTH/WIDTH/Length (mm)	Specification	Material
AA0081101	1.0/3/90	27 Holes	Pure Titanium
AA0081102	1.0/3/91	20Holes	

Note: Used with 1.1 cortical screw

» LOCKING PLATE SYSTEM

1.1mm Double T Cuttable Plate AA00911



Product Code	HEADWIDTH/WIDTH/Length (mm)	Specification	Material
AA0091101	1.0/3/90	30 Holes	Pure Titanium

Note: Used with 1.1mm cortical screw

1.5mm Double T Cuttable Plate AA01015



Product Code	HEADWIDTH/WIDTH/Length (mm)	Specification	Material
AA0101501	1.5/4/148	30 Holes	Pure Titanium
AA0101502	1.5/4/148	32 Holes	
AA0101503	1.5/4/152	32 Holes	

Note: Used with 1.5mm cortical screw

2.0mm Double T Cuttable Plate AA01120



Product Code	HEADWIDTH/WIDTH/Length (mm)	Specification	Material
AA0112001	1.8/5/160	29 Holes	Pure Titanium
AA0112002	1.8/5/160	31 Holes	
AA0112003	1.8/5/164	31 Holes	

Note: Used with 2.0mm cortical screw

2.4mm Double T Cuttable Plate AA01224



Product Code	HEADWIDTH/WIDTH/Length (mm)	Specification	Material
AA0122401	1.8/5/160	27 Holes	Pure Titanium
AA0122402	1.8/5/160	29 Holes	
AA0122403	1.8/5/164	29 Holes	

Note: Used with 2.4mm cortical screw

LOCKING PLATE SYSTEM

1.6mm Locking Screw AA00216



Product Code	HEADWIDTH/WIDTH/Length (mm)	Material
AA0021604-AA0021610	1.6x4/5/6/7/8/10	Titanium Alloy
Note: Match cross head screw driver, compatible for cortical and cancellous bone		

2.4mm Locking Screw AA00224



Product Code	HEADWIDTH/WIDTH/Length (mm)	Material
AA0022405-AA0022416	2.4x5/6/7/8/10/12/14/16	Titanium Alloy
Note: Match T6 torx screw driver, compatible for cortical and cancellous bone		

3.2mm Locking Screw AA00232



Product Code	HEADWIDTH/WIDTH/Length (mm)	Material
AA0023206-AA0023220	3.2x6/8/10/12/14/16/18/20	Titanium Alloy
Note: Match T6 torx screw driver, compatible for cortical and cancellous bone		

4.0mm Locking Screw AA00240



Product Code	HEADWIDTH/WIDTH/Length (mm)	Material
AA0024010-AA0024026	4.0x10/12/14/16/18/20/22/24/26	Titanium Alloy
Note: Match T10 torx screw driver, compatible for cortical and cancellous bone		

1.1mm Cortical Self-Tapping Screw AA00311



Product Code	HEADWIDTH/WIDTH/Length (mm)	Material
AA0031104-AA0031112	1.1x4/5/6/7/8/9/10/12	Titanium Alloy
Note: Match cross head screw driver, compatible for cortical and cancellous bone		

1.5mm Cortical Self-Tapping Screw AA00315



Product Code	HEADWIDTH/WIDTH/Length (mm)	Material
AA0031506-AA0031526	1.5x6/7/8/10/12/14/16/18/20/22/24/26	Titanium Alloy
Note: Match T6 torx screw driver, compatible for cortical and cancellous bone		

2.4mm Cortical Self-Tapping Screw AA00324



Product Code	HEADWIDTH/WIDTH/Length (mm)	Material
AA0032410-AA0032432	2.4x10/12/14/16/18/20/22/24/26/28/30/32	Titanium Alloy
Note: Match T10 torx screw driver, compatible for cortical and cancellous bone		


2.7mm Cortical Self-Tapping Screw AA00327



Product Code	HEADWIDTH/WIDTH/Length (mm)	Material
AA0032712-AA0032734	2.7x12/14/16/18/20/22/24/26/28/30/32/34	Titanium Alloy
Note: Match T10 torx screw driver, compatible for cortical and cancellous bone		


LOCKING PLATE SYSTEM

1.6mm Locking Plate Instrument Kit(A001)



NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION	QTY
1	A001-01	Drill Bit	Ø0.8mm	1
2	A001-02	Drill Bit	Ø1.1mm	1
3	A001-03	Locking Drill Guide	Ø1.1	1
4	A001-04	Drill Guide	Ø0.8/1.1mm	1
5	A001-05	Cross Head Screw Driver		1
6	A001-06	1.1mm T Type Plate(A Series)	2+3/2+4/2+5Holes	3pc;1pc/size
7	A001-07	1.1mm T Type Plate(B Series)	2+3/2+4/2+5Holes	3pc;1pc/size
8	A001-08	1.1mm T Type Plate(C Series)	3+3/3+4Holes	2pc;1pc/size
9	A001-09	1.1mm Straight Cuttable Plate	30Holes,90mm	1
10	A001-10	1.1mm Straight Cuttable Plate	23Holes,91mm	1
11	A001-11	1.1mm Resconstruction Cuttable Plate	30Holes,90mm	1
12	A001-12	1.1mm Neutral Cuttable Plate	27Holes,90mm	1
13	A001-13	1.6 Locking Plate	30Holes	1
14	A001-14	1.6 Locking Plate	26Holes	1
15	A001-15	Locking Screw	1.6x4/5/6/7/8/10mm	6pc/size
16	A001-16	Cortical Self Tapping Screw	1.1x4/5/6/7/8/9/10/12mm	6pc/size
17	A001-17	Implants Case		1


2.4mm Locking Plate Instrument Kit(A002)



NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION	QTY
1	A002-01	Drill Bit (AO)	Ø1.1mm	1
2	A002-02	Drill Bit (AO)	Ø1.5mm	1
3	A002-03	Drill Bit (AO)	Ø1.8mm	1
4	A002-04	Drill Bit Stopper	Ø1.8mm	1
5	A002-05	Drill Bit Stopper Spanner		1
6	A002-06	Locking Drill Guide	Ø1.1mm	1
7	A002-07	Locking Drill Guide	Ø1.8mm	1
8	A002-08	Compression Drill Guide	Ø1.1mm	1
9	A002-09	Drill Guide	Ø1.1/ Ø1.5mm	1
10	A002-10	AQ Quick Coupling Handle		1
11	A002-11	Depth Gauge	30mm	1
12	A002-12	Screw Holder		1
13	A002-13	AQ Quick Coupling Torx Handle Screw Driver	T6	1
14	A002-14	Plate Bending Wrench		1
15	A002-15	Implants Box		1
16	A002-16	Tool Box		1


LOCKING PLATE SYSTEM

3.2mm Locking Plate Instrument Kit(A003)



NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION	QTY
1	A003-01	Drill Bit (AO)	Ø1.8mm	1
2	A003-02	Drill Bit (AO)	Ø2.5mm	1
3	A003-03	Drill Bit Stopper	Ø2.5mm	1
4	A003-04	Drill Bit Stopper Spanner		1
5	A003-05	Locking Drill Guide	Ø1.8mm	1
6	A003-06	Compression Drill Guide	Ø1.8mm	1
7	A003-07	Locking Drill Guide	Ø2.5mm	1
8	A003-08	Drill Guide	Ø1.8/Ø2.5mm	1
9	A003-09	AQ Quick Coupling Handle		1
10	A003-10	Depth Gauge	50mm	1
11	A003-11	Screw Holder		1
12	A003-12	Torx Handle Screw Driver	T6	1
13	A003-13	Plate Bending Wrench		1
14	A003-14	Implants Box		1
15	A003-15	Tool Box		1

4.0mm Locking Plate Instrument Kit(A004)



NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION	QTY
1	A004-01	Drill Bit (AO)	Ø2.0mm	1
2	A004-02	Drill Bit (AO)	Ø2.7mm	1
3	A004-03	Drill Bit (AO)	Ø3.2mm	1
4	A004-04	Drill Bit Stopper	Ø3.2mm	1
5	A004-05	Drill Bit Stopper Spanner		1
6	A004-06	Locking Drill Guide	Ø2.0mm	1
7	A004-07	Locking Drill Guide	Ø3.2mm	1
8	A004-08	Compression Drill Guide	Ø2.0mm	1
9	A004-09	Drill Guide	Ø2.0/ Ø2.7mm	1
10	A004-10	AQ Quick Coupling Handle		1
11	A004-11	Depth Gauge	60mm	1
12	A004-12	Screw Holder		1
13	A004-13	AQ Quick Coupling Torx Handle Screw Driver	T10	1
14	A004-14	Plate Bending Wrench		1
15	A004-15	Implants Box		1
16	A004-16	Tool Box		1

LOCKING PLATE SYSTEM

2.0mm Straight Locking Plate AB00120



Product Code	Length (mm)	Specification	Material
AB0012004	23.2	4Holes	Pure Titanium
AB0012005	28.5	5Holes	
AB0012006	33.8	6Holes	
AB0012007	39.1	7Holes	
AB0012008	44.4	8Holes	
AB0012009	49.7	9Holes	
AB0012010	55	10Holes	

Note: Select Ø2.0 Locking Screw(Triangle) Ø2.0 Cortical Self Tapping Screw(Triangle), photo show 8 holes

2.0mm L Type Straight Locking Plate (L/R) AB00220



Product Code	Length (mm)	Specification	Material
AB0022005L/R	22.5	5 Holes L/R	Pure Titanium
AB0022006L/R	27.8	6 Holes L/R	
AB0022007L/R	33.1	7 Holes L/R	
AB0022008L/R	38.4	8 Holes L/R	

Note: Select Ø2.0 Locking Screw(Triangle) Ø2.0 Cortical Self Tapping Screw(Triangle), photo show 8 holes

2.0mm Condyle Locking Plate AB00320



Product Code	Length (mm)	Specification	Material
AB0032005	23.1	5 Holes	Pure Titanium
AB0032006	28.4	6 Holes	
AB0032007	33.7	7 Holes	
AB0032008	39	8 Holes	
AB0032009	44.3	9 Holes	

Note: Select Ø2.0 Locking Screw(Triangle) Ø2.0 Cortical Self Tapping Screw(Triangle), photo show 8 holes

2.0mm T Type Locking Plate AB00420



Product Code	Length (mm)	Specification	Material
AB0042005	23.4	5 Holes	Pure Titanium
AB0042006	28.7	6 Holes	
AB0042007	34	7 Holes	
AB0042008	39.3	8 Holes	
AB0042009	44.6	9 Holes	

Note: Select Ø2.0 Locking Screw(Triangle) Ø2.0 Cortical Self Tapping Screw(Triangle), photo show 8 holes

2.5mm Straight Locking Plate AB00525



Product Code	Length (mm)	Specification	Material
AB0052504	30	4 Holes	Pure Titanium
AB0052505	37	5 Holes	
AB0052506	44	6 Holes	
AB0052507	51	7 Holes	
AB0052508	58	8 Holes	
AB0052509	65	9 Holes	
AB0052510	72	10 Holes	

Note: Select Ø2.5 Locking Screw Ø2.5 Cortical Self Tapping Screw, photo show 10 holes

» LOCKING PLATE SYSTEM

2.5mm L Type Locking Plate(L/R) AB00625



Product Code	Length (mm)	Specification	Material
AB0062505L/R	29.5	5Holes L/R	Pure Titanium
AB0062506L/R	36.5	6Holes L/R	
AB0062507L/R	43.5	7Holes L/R	
AB0062508L/R	50.5	8Holes L/R	
AB0062509L/R	57.5	9Holes L/R	

Note: Select Ø2.5 Locking Screw, Ø2.5 Cortical Self Tapping Screw, photo show 9 holes

2.5mm Condyle Locking Plate AB00725



Product Code	Length (mm)	Specification	Material
AB0072505	30.5	5Holes	Pure Titanium
AB0072506	37.5	6Holes	
AB0072507	44.5	7Holes	
AB0072508	51.5	8Holes	
AB0072509	58.5	9Holes	

Note: Select Ø2.5 Locking Screw, Ø2.5 Cortical Self Tapping Screw, photo show 9 holes

2.5mm L Type Locking Plate AB00825



Product Code	Length (mm)	Specification	Material
AB0082505	30.5	5Holes	Pure Titanium
AB0082506	37.5	6Holes	
AB0082507	44.5	7Holes	
AB0082508	51.5	8Holes	
AB0082509	58.5	9Holes	

Note: Select Ø2.5 Locking Screw, Ø2.5 Cortical Self Tapping Screw, photo show 9 holes

2.7mm Locking Plate AB00927



Product Code	Length (mm)	Specification	Material
AB0092704	54	4Holes	Pure Titanium
AB0092705	64	5Holes	
AB0092706	74	6Holes	
AB0092707	84	7Holes	
AB0092708	94	8Holes	
AB0092709	104	9Holes	
AB0092710	114	10Holes	

Note: Select Ø2.7 Locking Screw, Ø2.7 Cortical Self Tapping Screw, photo show 7 holes

2.7mm Reconstruction Locking Plate AB01027



Product Code	Length (mm)	Specification	Material
AB0102704	40.5	4Holes	Pure Titanium
AB0102705	50.5	5Holes	
AB0102706	60.5	6Holes	
AB0102707	70.5	7Holes	
AB0102708	80.5	8Holes	
AB0102709	90.5	9Holes	
AB0102710	100.5	10Holes	

Note: Select Ø2.7 Locking Screw, Ø2.7 Cortical Self Tapping Screw, photo show 7 holes

» LOCKING PLATE SYSTEM

2.0, 2.5, 2.7mm Locking Screw AC00120, AC00125, AC00127



Product Code	Diameter(mm)	Length(mm)	Material
AC0012006-AC0012022	29.5	5Holes L/R	Titanium Alloy
AC0012506-AC0012530	36.5	6Holes L/R	
AC0012706-AC0012730	43.5	7Holes L/R	

Note: 2.0 Locking screw match triangle screw driver, 2.5/2.7 Locking screw match T2.3 screw driver


2.0, 2.5, 2.7mm Cortical Self Tapping Screw AD00120, AD00125, AD00127



Product Code	Diameter(mm)	Length(mm)	Material
AD0012006-AC0012022	2.0(Triangle)	6-22 2mm size gap	Titanium Alloy
AD0012506-AC0012530	2.5(Torx)	6-30 2mm size gap	
AD0012706-AC0012730	2.7(Torx)	6-30 2mm size gap	

Note: 2.0 cortical self tapping screw match triangle screw driver, 2.5/2.7 cortical self tapping screw match T2.3 screw driver

2.0 Locking Plate Instrument Kit(A005)



NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION
1	A005-01	Triangle Screw Driver	
2	A005-02	Locking Screw Tapper	HA2
3	A005-03	Cortical Self Tapping Tapper	HA2
4	A005-04	Drill Bit	Ø1.5mm
5	A005-05	Drill Guide	Ø1.5mm
6	A005-06	Depth Gauge	30mm
7	A005-07	Screw Box	

2.5 Locking Plate Instrument Kit(A006)



NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION
1	A006-01	Torx Head Screw Driver	T2.3
2	A006-02	Locking Screw Tapper	HA2.5
3	A006-03	Cortical Self Tapping Tapper	HA2.5
4	A006-04	Drill Bit	Ø2.0mm
5	A006-05	Drill Guide	Ø2.0mm
6	A006-06	Depth Gauge	30mm
7	A006-07	Screw Box	

2.7 Locking Plate Instrument Kit(A007)



NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION
1	A007-01	Torx Head Screw Driver	T2.3
2	A007-02	Locking Screw Tapper	HA2.7
3	A007-03	Cortical Self Tapping Tapper	HA2.7
4	A007-04	Drill Bit	Ø2.1mm
5	A007-05	Drill Guide	Ø2.1mm
6	A007-06	Depth Gauge	60mm
7	A007-07	Screw Box	

NON LOCKING PLATE SYSTEM

1.5mm Straight Plate AF0011520



Product Code	Thickness/Width Length (mm)	Specification	Material
AF0011520	1.2/3.8/100	20Holes	Stainless Steel

Note: Photo for 20 holes, select Ø1.5mm cortical self tapping screw

1.5mm Straight Plate(Thickening) AF00215030519



Product Code	Thickness/Width Length (mm)	Specification	Material
AF00215030519	1.5/3.8/150	3+5+19Holes	Stainless Steel

Note: Photo for 3+5+19 holes, select Ø1.5mm cortical self tapping screw

2.0mm Straight Compression Plate Type I AF00320



Product Code	Thickness/Width Length (mm)	Specification	Material
AF0032006	1.5/5/33	6Holes	Stainless Steel
AF0032007	1.5/5/38	7Holes	
AF0032008	1.5/5/43	8Holes	
AF0032009	1.5/5/48	9Holes	
AF0032010	1.5/5/53	10Holes	
AF0032012	1.5/5/63	12Holes	

Note: Photo for 10 holes, select Ø2.0mm cortical self tapping screw

2.0mm Straight Compression Plate Type II AF00420



Product Code	Thickness/Width Length (mm)	Specification	Material
AF0042004	1.2/5.0/23	4Holes	Stainless Steel
AF0042005	1.2/5.0/28	5Holes	
AF0042006	1.2/5.0/33	6Holes	

Note: Photo for 6 holes, select Ø2.0mm cortical self tapping screw

2.0mm Straight Round Hole Plate Type I AF00520



Product Code	Thickness/Width Length (mm)	Specification	Material
AF0052004	1.2/5/26	4Holes	Stainless Steel
AF0052006	1.2/5/38	6Holes	
AF0052008	1.2/5/50	8Holes	
AF0052010	1.2/5/62	10Holes	

Note: Photo for 10 holes, select Ø2.0mm cortical self tapping screw

NON LOCKING PLATE SYSTEM

2.0mm Straight Round Hole Plate Type II AF00620



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0062004	1.2/5/31	4Holes	Stainless Steel
AF0062005	1.2/5/36	5Holes	

Note: Photo for 5 holes, select Ø2.0mm cortical self tapping screw

2.0mm Straight Round Hole Plate Type III AF00720



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0072004	1.2/5/26	4Holes	Stainless Steel
AF0072005	1.2/5/31	5Holes	

Note: Photo for 5 holes, select Ø2.0mm cortical self tapping screw

2.0mm Straight Round Hole Plate Type VI AF00820



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0082006	1.5/5.5/50	6Holes	Stainless Steel
AF0082007	1.5/5.5/60	7Holes	
AF0082008	1.5/5.5/70	8Holes	

Note: Photo for 7 holes, select Ø2.0mm cortical self tapping screw

2.0mm Straight Plate AF00920



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0092006	1.5/5/41	6Holes	Stainless Steel
AF0092008	1.5/5/52	8Holes	

Note: Photo for 8 holes, select Ø2.0mm cortical self tapping screw

2.0mm Straight Neutral Plate AF01020



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0102006	1.5/5/40	6Holes	Stainless Steel
AF0102008	1.5/5/52	8Holes	
AF0102010	1.5/5/64	10Holes	

Note: Photo for 8 holes, select Ø2.0mm cortical self tapping screw

» NON LOCKING PLATE SYSTEM

2.0mm Straight Plate(Cuttable) AF01120



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0112050	1.2/5/250	50Holes	Stainless Steel

Note: Photo for 50 holes, select Ø2.0mm cortical self tapping screw

2.4mm LC-DCP AF01224



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0122404	2.0/6.5/34	4Holes	Stainless Steel
AF0122405	2.0/6.5/42	5Holes	
AF0122406	2.0/6.5/50	6Holes	
AF0122407	2.0/6.5/58	7Holes	
AF0122408	2.0/6.5/66	8Holes	
AF0122410	2.0/6.5/82	10Holes	

Note: Photo for 8 holes, select Ø2.4mm cortical self tapping screw

2.4mm Straight Neutral Plate AF01324



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0132406	2.0/6.5/52.5	6Holes	Stainless Steel
AF0132408	2.0/6.5/65.5	8Holes	

Note: Photo for 6 holes, select Ø2.4mm cortical self tapping screw

2.4mm Straight Neutral/Lenghtening Plate AF01424



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0142406	2.0/6.5/59	6Holes	Stainless Steel
AF0142408	2.0/6.5/72	8Holes	

Note: Photo for 8 holes, select Ø2.4mm cortical self tapping screw

2.4mm Straight Plate(Cuttable) AF01524



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0152422	1.5/6.5/162.5	22Holes	Stainless Steel

Note: Photo for 22 holes, select Ø2.4mm cortical self tapping screw

» NON LOCKING PLATE SYSTEM

2.7mm Straight Compression Plate Type I AF01627



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0162704	2.0/8/37	4Holes	Stainless Steel
AF0162705	2.0/8/45	5Holes	
AF0162706	2.0/8/53	6Holes	
AF0162707	2.0/8/61	7Holes	
AF0162708	2.0/8/69	8Holes	
AF0162709	2.0/8/77	9Holes	
AF0162710	2.0/8/85	10Holes	
AF0162712	2.0/8/101	12Holes	

Note: Photo for 7 holes, select Ø2.7mm cortical self tapping screw

2.7mm LC-DCP AF01927



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0192704	2.0/8/43	4Holes	Stainless Steel
AF0192705	2.0/8/52	5Holes	
AF0192706	2.0/8/61	6Holes	
AF0192707	2.0/8/70	7Holes	
AF0192708	2.0/8/79	8Holes	
AF0192710	2.0/8/98	10Holes	
AF0192712	2.0/8/117	12Holes	

Note: Photo for 10 holes, select Ø2.7mm cortical self tapping screw

2.7mm Round Hole Neutral Plate AF02027



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0202708	2.5/8.0/97	8Holes	Stainless Steel

Note: Photo for 8 holes, select Ø2.7mm cortical self tapping screw

3.5mm Straight Compression Plate AF02135



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0213506	3.0/10/77	6Holes	Stainless Steel
AF0213508	3.0/10/101	8Holes	
AF0213510	3.0/10/125	10Holes	
AF0213512	3.0/10/149	12Holes	

Note: Photo for 8 holes, select Ø3.5mm cortical self tapping screw

» NON LOCKING PLATE SYSTEM

3.5mm Straight Compression Plate(Wide) AF02235



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0223506	3.5/12/88	6Holes	Stainless Steel
AF0223508	3.5/12/116	8Holes	
AF0223510	3.5/12/144	10Holes	
AF0223512	3.5/12/172	12Holes	
AF0223514	3.5/12/200	14Holes	

Note: Photo for 8 holes, select Ø3.5mm cortical self tapping screw

1.5mm Straight Reconstruction Plate AF0231520



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0231520	1.3/4/100	20Holes	Stainless Steel

Note: Photo for 20 holes, select Ø1.5mm cortical self tapping screw

2.0mm Straight Reconstruction Plate AF0242020



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0242020	1.5/5/100	20Holes	Stainless Steel

Note: Photo for 20 holes, select Ø2.0mm cortical self tapping screw

2.4mm Straight Reconstruction Plate AF0252425



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0252425	1.5/6.5/162.5	25Holes	Stainless Steel

Note: Photo for 25 holes, select Ø2.4mm cortical self tapping screw

2.7mm Straight Reconstruction Plate AF0262720



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0262720	2.5/8/160	20Holes	Stainless Steel

Note: Photo for 20 holes, select Ø2.7mm cortical self tapping screw

» NON LOCKING PLATE SYSTEM

3.5mm Straight Reconstruction Plate AF0273520



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0273520	3.0/10/144	12Holes	Stainless Steel

Note: Photo for 12 holes, select Ø3.5mm cortical self tapping screw

1.5mm Double T Type Plate AF0281530



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0281530	1.5/4/154	30Holes	Stainless Steel

Note: Photo for 30 holes, select Ø1.5mm cortical self tapping screw

1.5mm T Type Plate(3 Holes in Head) AF029150304



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF029150304	1.2/4/29	3+4Holes	Stainless Steel

Note: Photo for 3+4 holes, select Ø1.5mm cortical self tapping screw

2.0mm T Type Plate(3 Holes in Head) AF030200308



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF030200308	1.5/5.0/62	3+8Holes	Stainless Steel

Note: Photo for 3+8 holes, select Ø2.0mm cortical self tapping screw

2.0mm T Type Plate(2 Holes in Head)-I AF03120



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF031200204	1.2/5.0/38	2+4Holes	Stainless Steel
AF031200208	1.2/5.0/56.5	2+8Holes	

Note: Photo for 2+4 holes, select Ø2.0mm cortical self tapping screw

» NON LOCKING PLATE SYSTEM

2.0mm T Type Plate(2 Holes in Head)-II AF032200204



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF032200204	1.5/5/42	2+4Holes	Stainless Steel

Note: Photo for 2+4 holes, select Ø2.0mm cortical self tapping screw

2.4mm T Type Plate(2 Holes in Head) AF033240208



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF033240208	1.5/6.5/60	2+8Holes	Stainless Steel

Note: Photo for 2+8 holes, select Ø2.4mm cortical self tapping screw

2.7mm T Type Plate(2 Holes in Head) AF03427



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF034270204	2.0/8/41	2+4Holes	Stainless Steel
AF034270208	2.0/8/97.5	2+8Holes	

Note: Photo for 2+8 holes, select Ø2.7mm cortical self tapping screw

3.5mm T Type Plate(2 Holes in Head) AF035350208



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF035350208	3.0/11/97.5	2+8Holes	Stainless Steel

Note: Photo for 2+8 holes, select Ø3.5mm cortical self tapping screw

2.0mm T Type Round Hole Plate AF03620



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF036200204L	1.2/5/21/90°	2+4Holes L	Stainless Steel
AF036200204R		2+4Holes R	

Note: Photo for 2+4 holes, select Ø2.0mm cortical self tapping screw

2.0mm 120° L Type Round Hole Plate AF03720



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF037200204L	1.2/5/24.5/120°	2+4Holes L	Stainless Steel
AF037200204R		2+4Holes R	

Note: Photo for 2+4 holes, select Ø2.0mm cortical self tapping screw

» NON LOCKING PLATE SYSTEM

2.0mm L Type Curved Plate AF03820



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF03820L	1.5/10/26	2+4Holes L	Stainless Steel
AF03820R		2+4Holes R	

Note: Photo for 2+4 holes, select Ø2.0mm cortical self tapping screw

2.7mm L Type Curved Plate AF03927



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF03927L	2.0/16.3/40.7	2+4Holes L	Stainless Steel
AF03927R		2+4Holes R	

Note: Photo for 2+4 holes, select Ø2.7mm cortical self tapping screw

3.5mm L Type Curved Plate AF04035



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF04035L	2.5/16.9/63.8	2+4Holes L	Stainless Steel
AF04035R		2+4Holes R	

Note: Photo for 2+4 holes, select Ø3.5mm cortical self tapping screw

2.0mm Distal Femoral Compression Plate AF04120



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF04120L	1.5/5.5/75.5	16Holes L	Stainless Steel
AF04120R		16Holes R	

Note: Photo for 16 holes, select Ø2.0mm cortical self tapping screw

2.0mm Distal Femoral Round Hole Plate AF04220



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0422008L	1.5/5.5/43.7	8Holes L	Stainless Steel
AF0422008R		8Holes R	

Note: Photo for 8 holes, select Ø2.0mm cortical self tapping screw

2.7mm Distal Femoral Plate Type I AF04327



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF04327L	2.5/8/68	3+7Holes L	Stainless Steel
AF04327R		3+7Holes R	

Note: Photo for 3+7 holes, select Ø2.7mm cortical self tapping screw

» NON LOCKING PLATE SYSTEM

2.7mm Distal Femoral Plate Type II AF04427



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF04427L	2.5/8/68	4+5Holes L	Stainless Steel
AF04427R		4+5Holes R	
Note: Photo for 4+5 holes, select Ø2.7mm cortical self tapping screw			

3.5mm Distal Femoral Plate AF04535



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF04535L	3.0/11/115	3+10Holes L	Stainless Steel
AF04535R		3+10Holes R	
Note: Photo for 3+10 holes, select Ø3.5mm cortical self tapping screw			

1.5mm Acetabular Plate AF0461504



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0461504	1.2/4/19.1	4Holes	Stainless Steel
Note: Photo for 4 holes, select Ø1.5mm cortical self tapping screw			

2.0mm Acetabular Plate AF0472004



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0472004	1.2/5/26.8	4Holes	Stainless Steel
Note: Photo for 4 holes, select Ø2.0mm cortical self tapping screw			

2.4mm Acetabular Plate AF0482406



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0482406	1.5/5/33.4	6Holes	Stainless Steel
Note: Photo for 6 holes, select Ø2.4mm cortical self tapping screw			

2.7mm Acetabular Plate AF0492706



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0492706	2.0/7/40	6Holes	Stainless Steel
Note: Photo for 6 holes, select Ø2.7mm cortical self tapping screw			

» NON LOCKING PLATE SYSTEM

3.5mm Acetabular Plate AF0503506



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0503506	2.5/10/48	6Holes	Stainless Steel
Note: Photo for 6 holes, select Ø3.5mm cortical self tapping screw			

1.5mm X Type Plate(Cuttable) AF05115



Product Code	Thickness/Width/Length(mm)	Material
AF0511501	1.2/9/14	Stainless Steel
AF0511502	1.2/11/15	
AF0511503	1.2/11/15	
Note: Photo for 4 holes, select Ø1.5mm cortical self tapping screw		

2.0mm X Type Plate(Cuttable) AF05220



Product Code	Thickness/Width/Length(mm)	Material
AF0522001	1.2/10.4/15.8	Stainless Steel
AF0522002	1.2/12.9/17	
AF0522003	1.2/10.6/22.9	
Note: Photo for 4 holes, select Ø2.0mm cortical self tapping screw		

2.7mm X Type Plate(Cuttable) AF05327



Product Code	Thickness/Width/Length(mm)	Material
AF0532701	2.0/13.9/20.9	Stainless Steel
AF0532702	2.0/16.4/24.8	
AF0532703	2.0/13.9/29.8	
Note: Photo for 6 holes, select Ø2.7mm cortical self tapping screw		

Model For Cutting



1.5/2.0mm Fusion Plate AF0541513



Product Code	Thickness/Width/Length(mm)	Specification	Material
AF0541513	1.3/6/80	13Holes	Stainless Steel
Note: Photo for 13 holes, select Ø2.0mm cortical self tapping screw			

NON LOCKING PLATE SYSTEM

Ø1.5, 2.0, 2.4, 2.7, 3.5mm Cortical Self-Tapping Screw(Hex)
AG00115, AG00220, AG00324, AG00427, AG00535



Product Code	Diameter(mm)	Length(mm)	Material
AG0011506-AG0011522	HA1.5	6-22 Every 2mm is a specification	Stainless Steel
AG0022006-AG0022022	HA2.0	6-22 Every 2mm is a specification	
AG0032406-AG0032426	HA2.4	6-26 Every 2mm is a specification	
AG0042706-AG0042730	HA2.7	8-30 Every 2mm is a specification	
AG0053506-AG0053536	HA3.5	6-36 Every 2mm is a specification	

Note: HA1.5, HA2.0, HA2.4 select SW1.5 Hex screw driver, HA2.7, HA3.5 select SW2.5 Hex screw driver

Screw Sterilization Box 1.5/2.0/2.7 (A008)



Product Code	No.	Description	Specification	Quantities
A008-01	1	1.5 Self-Tapping Screw(Hex)	1.5x6/8/10/12/14/16/18/20mm	5pcs for each size
A008-02	2	2.0 Self-Tapping Screw(Hex)	2.0x8/10/12/14/16/18/20/22/24/26/28mm	6pcs for each size
A008-03	3	2.7 Self-Tapping Screw(Hex)	2.7x10/12/14/16/18/20/22/24/26/28/30mm	6pcs for each size





Screw Sterilization Box 1.5/2.0/2.4/2.7/3.5 (A009)



Product Code	No.	Description	Specification	Quantities
A009-01	1	1.5 Self-Tapping Screw(Hex)	1.5x6/8/10/12/14/16/18/20mm	5pcs for each size
A009-02	2	2.0 Self-Tapping Screw(-Hex)	2.0x8/10/12/14/16/18/20/22/24/26/28mm	5pcs for each size
A009-03	3	2.4 Self-Tapping Screw(-Hex)	2.4x8/10/12/14/16/18/20/22/24/26/28/30mm	5pcs for each size
A009-04	4	2.7 Self-Tapping Screw(-Hex)	2.7x10/12/14/16/18/20/22/24/26/28/30mm	5pcs for each size
A009-05	5	3.5 Self-Tapping Screw(-Hex)	3.5x8/10/12/14/16/18/20/22/24/26/28/30mm	5pcs for each size

NON LOCKING PLATE SYSTEM

Veterinary Plate Set (A010)



NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION	NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION
1	AF0262720	3.5mm Straight Reconstruction Plate	12 holes	25	A010-25	2.4mm Distal Femoral Compression Plate	12 holes L
2	AF0262720	2.7mm Straight Reconstruction Plate	20 holes	26	A010-26	2.4mm Distal Femoral Compression Plate	12 holes R
3	AF0252425	2.4mm Straight Reconstruction Plate	25 holes	27	AF04327L	2.7mm Distal Femoral Compression Plate	3+7 holes L
4	AF0242020	2.0mm Straight Reconstruction Plate	20 holes	28	AF04327R	2.7mm Distal Femoral Compression Plate	3+7 holes R
5	AF0231520	1.5mm Straight Reconstruction Plate	20 holes	29	AF04535L	3.5mm Distal Femoral Compression Plate	3+10 holes L
6	AF0032010	2.0mm Straight Compression Plate Type-I	10 holes	30	AF04535R	3.5mm Distal Femoral Compression Plate	3+10 holes R
7	AF0032008	2.0mm Straight Compression Plate Type-I	8 holes	31	A010-31	2.0mm DCP	9 holes
8	AF0122408	2.4mm LC-DCP	8 holes	32	A010-32	2.0mm DCP	8 holes
9	AF00122410	2.4mm LC-DCP	10 holes	33	A010-33	2.4mm DCP	12 holes
10	AF0162710	2.7mm Straight Compression Plate Type-I	10 holes	34	A010-34	2.4mm DCP	14 holes
11	AF0162712	2.7mm Straight Compression Plate Type-I	12 holes	35	A010-35	2.7mm DCP	9 holes
12	AF0461504	1.5mm Acetabular Plate	4 holes	36	A010-36	2.7mm DCP	8 holes
13	AF0472004	2.0mm Acetabular Plate	4 holes	37	A010-37	3.5mm DCP	12 holes
14	AF0482406	2.4mm Acetabular Plate	6 holes	38	A010-38	3.5mm DCP	14 holes
15	AF0492706	2.7mm Acetabular Plate	6 holes	39	A010-39	1.5mm Cuttable Plate	27 holes
16	AF0503506	3.5mm Acetabular Plate	6 holes	40	A010-40	1.5mm Double T Type Plate	30 holes
17	AF029150304	1.5mm T Type Plate	3+4 holes	41	A010-41	2.0mm Cuttable Plate	50 holes
18	AF032200204	2.0mm T Type Plate-II	2+4 holes	42	A010-42	2.4mm Cuttable Plate	22 holes
19	AF033240208	2.4mm T Type Plate	2+8 holes	43	A010-43	2.0mm Neutral Plate	6 holes
20	AF0142408	2.4mm Straight Neutral Plate	8 holes	44	A010-44	2.0mm Neutral Plate	7 holes
21	AF0102008	2.0mm Straight Neutral Plate	8 holes	45	A010-45	2.4mm Neutral Plate	6 holes
22	AF0202708	2.7mm Round Holle Neutral Plate	8 holes	46	A010-46	2.7mm T Type Plate	2+8 holes
23	AF04120L	2.0mm Distal Femoral Compression Plate	16 holes L	47	A010-47	3.5mm T Type Plate	2+8 holes
24	AF04120R	2.0mm Distal Femoral Compression Plate	16 holes R	48	A010-48	Bone Plate Box	

Sterilization Box A011






Product Code: A011
Specification: 275x165x50mm

Sterilization Box A012

Product Code: A012
Specification: 300x180x105mm

NON LOCKING PLATE SYSTEM

Veterinary Orthopedic Instrument 1.5/2.0/2.4/2.7/3.5mm (A013)



NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION
1	A013-01	Reduction Forceps with point-Medium	160mm
2	A013-02	Reduction Forceps with point-Small	135mm
3	A013-03	Self-centering bone holding Forcep	X-Small
4	A013-04	Reduction Forcep	180mm
5	A013-05	Self-centering bone holding Forcep	Small
6	A013-06	Wire Cutter	2.5mm
7	A013-07	End Cutting Forcep	2.5mm
8	A013-08	Wire Guider	Ø20mm
9	A013-09	K Wire punch	1.5-2.5mm
10	A013-10	K Wire Bender	
11	A013-11	Ruskin Ronger	2mm
12	A013-12	Bone Rasp	
13	A013-13	Double End Periosteal Elevator	5mm
14	A013-14	Mini Hohman Retractor	6mm
15	A013-15	Mini Hohman Retractor	8mm
16	A013-16	Mini Hohman Retractor	12mm
17	A013-17	Screw Holder	Small
18	A013-18	AO Plates Bender	1.5/2.0/2.4mm
19	A013-19	AO Plates Bender	2.7/3.5mm
20	A013-20	Depth Gauge	0-30mm
21	A013-21	Depth Gauge	0-40mm
22	A013-22	Screw Driver	SW1.5
23	A013-23	Screw Driver	SW2.5
24	A013-24	AO Quick Couping T Type Handle	
25	A013-25	AO Quick Couping Countersink	Ø4mm

NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION
26	A013-26	AO Quick Couping Countersink	Ø6mm
27	A013-27	AO Quick Couping Screw Tap	1.5mm
28	A013-28	AO Quick Couping Screw Tap	2.0mm
29	A013-29	AO Quick Couping Screw Tap	2.4mm
30	A013-30	AO Quick Couping Screw Tap	2.7mm
31	A013-31	AO Quick Couping Screw Tap	3.5mm
32	A013-32	Double Drill Guider	1.1/1.5mm
33	A013-33	Double Drill Guider	1.5/2.0mm
34	A013-34	Double Drill Guider	1.8-2.4mm
35	A013-35	Double Drill Guider	2.0/2.7mm
36	A013-36	Double Drill Guider	2.5/3.5mm
37	A013-37	Drill Sleeve	1.5mm
38	A013-38	Drill Sleeve	1.8mm
39	A013-39	Drill Sleeve	2.0mm
40	A013-40	Drill Sleeve	2.5mm
41	A013-41	Bone Plate Model	4pc
42	A013-42	Skin Retractor	
43	A013-43	Self Retaning Retractor	140mm
44	A013-44	Drill Bit	1.1mm
45	A013-45	Drill Bit	1.5mm
46	A013-46	Drill Bit	1.8mm
47	A013-47	Drill Bit	2.0mm
48	A013-48	Drill Bit	2.5mm
49	A013-49	Drill Bit	2.7mm

NON LOCKING PLATE SYSTEM

Suture-Wire AH001




Product Code	Diameter(mm)	Material
AH00104	0.4	Stainless Steel
AH00106	0.6	
AH00108	0.8	
AH00110	1.0	
AH00112	1.2	

Kirschner Wire(Single/Double Head) AH002



Product Code	Diameter(mm)	Length	Material
AH00208	0.8	250	Stainless Steel
AH00210	1.0	250	
AH00212	1.2	250	
AH00215	1.5	250	
AH00220	2.0	250	
AH00225	2.5	250	
AH00230	3.0	250	
AH00235	3.5	250	
AH00240	4.0	250	
AH00245	4.5	250	

K Wire Racks AH003



Product Code	Specification(mm)	Description
AH003	330x140x50mm	used for 0.8/1.0/1.2/1.5/2.0/2.5/3.0/3.5/4.0/4.5 wire

Temporary Fixation Wire For Bone Plate AH004



Product Code	Diameter(mm)	Length	Material
AH00401	1.2	85	Stainless Steel
AH00402	1.5	150	

Note: 1.2mm used for 1.5/2.0/2.4mm bone plate 1.5mm used for 2.7/3.5mm bone plate

NON LOCKING PLATE SYSTEM

Drill Bit AZ001



Product Code	Diameter(mm)	Length(mm)	Material
AZ00108	0.8	60	Stainless Steel
AZ00111	1.1	85	
AZ00112	1.2	115	
AZ00115	1.5	115	
AZ00118	1.8	115	
AZ00120	2.0	115	
AZ00125	2.5	115	
AZ00127	2.7	115	
AZ00132	3.2	150	

AO Quick Coupling Drill Bit AZ002



Product Code	Diameter(mm)	Length(mm)	Material
AZ00208	0.8	60	Stainless Steel
AZ00211	1.1	85	
AZ00215	1.5	100	
AZ00218	1.8	125	
AZ00220	2.0	125	
AZ00225	2.5	150	
AZ00227	2.7	150	
AZ00232	3.2	150	

Elastic Nail AZ003



Product Code	Diameter(mm)	Length(mm)	Material
AZ00312	1.2	250	Titanium Alloy
AZ00315	1.5	250	
AZ00318	1.8	250	
AZ00320	2.0	250	
AZ00325	2.5	250	
AZ00330	3.0	250	




End Cap AZ004



Product Code	Specification	Description	Material
AZ004-01	small	used for 1.2/1.5/1.8/2.0mm Elastic Nail	Titanium Alloy
AZ004-02	big	used for 2.5/3.0mm Elastic Nail	


NON LOCKING PLATE SYSTEM

Veterinary Elastic Nail Instrument A0014




NO.	PRODUCT CODE	ITEM	SPECIFICATION
1	A0014-01	AO elastic nail	1.2x250mm
2	A0014-02	AO elastic nail	1.5x250mm
3	A0014-03	AO elastic nail	1.8x250mm
4	A0014-04	AO elastic nail	2.0x250mm
5	A0014-05	AO elastic nail	2.5x250mm
6	A0014-06	AO elastic nail	3.0x250mm
7	A0014-07	K wire bender	
8	A0014-08	K wire impactor	
9	A0014-09	AWL	
10	A0014-10	Cannulated hand drill	1-handle
11	A0014-11	Sliding hammer	
12	A0014-12	Wire cutter	Max 3.0mm, 230mm
13	A0014-13	Pin forcep	
14	A0014-14	Sliding hammer guide rod	
15	A0014-15	Sterilization box	

Veterinary Trochlea Rectangular Sulcoplasty Instruments Set A015




NO.	PRODUCT CODE	DESCRIPTION	SPECIFICATION
1	A015-01	Bone saw 0.4mm	140mm
2	A015-02	Small trochlea rasp	Ø3.5mm
3	A015-03	mini trochlea rasp	Ø2.2mm
4	A015-04	Bone chisel	3mm
5	A015-05	Bone chisel	4mm
6	A015-06	Bone chisel	5mm
7	A015-07	Bone chisel	6mm
8	A015-08	Bone chisel	8mm
9	A015-09	Bone chisel	10mm
10	A015-10	Bone chisel	15mm
11	A015-11	Double head bone rasp	
12	A015-12	Bone hammer	270g
13	A015-13	Sterilization box	



Aiming Reduction Forceps

Product Code	Specification
BA0010130	Small, Max 30mm, Length 144mm Small aiming device 1.1/1.5/1.8mm drill sleeve
BA0010260	Large, Max 60mm, Length 178mm Large aiming device 1.8/2.0/2.5mm drill sleeve



Multi-Functional Bone Holding Forceps

Product Code	Specification
BA0020130	Max 20mm, Length 139mm





Plate Holding Forceps

Product Code	Specification
BA0030134	Max 34mm, Length 162mm



Reduction Forceps With Point

Product Code	Specification
BA004140	Small, Max 30mm, Length 140mm
BA004176	Large, Max 62mm, Length 176mm



Reduction Forceps With Limit

Product Code	Specification
BA0050164	Max 64mm, Length 175mm



AO Self Locking Reduction Forcep

Product Code	Specification
BA006150	Length 150mm




AO Self Locking Reduction Forcep With Point

Product Code	Specification
BA007150	Length 150mm



Reduction Forcep With Twin Point

Product Code	Specification
BA008180	Length 180mm



Reduction Forceps With Point

Product Code	Specification
BA009135	Small Length 135mm
BA009160	Medium Length 160mm



AO Reduction Forceps With Jaws

Product Code	Specification
BA010140	Length 140mm



Reduction Forceps With Jaws

Product Code	Specification
BA011150	Length 150mm



Reduction Forcep

Product Code	Specification
BA012180	Length 180mm



BW Reduction Forcep

Product Code	Specification
BA013150	Length 150mm



BW Bone Holding Forcep

Product Code	Specification
BA014150	Length 150mm



MINI Bone Holding Forcep

Product Code	Specification
BA0150314	Holding 3mm, Length 140mm
BA0150514	Holding 5mm, Length 140mm
BA0150714	Holding 7mm, Length 140mm



Self-Centering Bone Holding Forcep

Product Code	Specification
BA0160115	Mini Length 150mm
BA0160218	Small Length 180mm



KERN Bone Holding Forcep

Product Code	Specification
BA0170114	Mini Length 140mm
BA0171218	Small Length 180mm



Lowman Bone Clamp Veterinary Pattern Cat Size-Small

Product Code	Specification
BA0180108	Small<15mm Length 180mm
BA0180212	Big<30mm Length 120mm





Plate Bone Forceps

Product Code	Specification
BA019160	Length 160mm



Pin Cutting Forcep

Product Code	Specification
BA020205	Max 2.5mm, Length 205mm



Wire Cutter

Product Code	Specification
BA0216053	Max 6mm, Length 530mm



Wire Cutter

Product Code	Specification
BA0224037	Max 4mm, Length 370mm



Front End Cutting Forcep

Product Code	Specification
BA0232523	Max 2.5mm, Length 235mm



Wire Plier And Cutter

Product Code	Specification
BA024165	Length 165mm



Wire Cutter

Product Code	Specification
BA0252521	Max 2.5mm, Length 210mm



Wire Pliers(Flat)

Product Code	Specification
BA026220	Length 220mm



Wire Pliers(Sharp)

Product Code	Specification
BA027200	Length 200mm



Wire/Small Pin Cutter

Product Code	Specification
BA028150	Max 1.8mm, Length 180mm



Wire/Small Pin Cutter

Product Code	Specification
BA0291214	Max 1.2mm, Length 140mm



Plate/Pin Cutting Plier TC Gold

Product Code	Specification
BA0302823	Max 2.8mm, Length 230mm



Plate/Pin Cutting Plier TC Gold

Product Code	Specification
BA0312217	Max 2.2mm, Length 170mm
BA0312823	Max 2.8mm, Length 230mm



Front End Wire Cutting Plier TC Gold

Product Code	Specification
BA0322823	Max 2.8mm, Length 230mm



Jacobs Chuck(Straight)

Product Code	Specification
BA03301	Straight Handle 0-4mm



Jacobs Chuck(T handle)

Product Code	Specification
BA03302	1-handle 0-6mm



Wire Guider

Product Code	Specification
BA03405	Small, Dia 13.8mm, Hole Dia, 0.5mm
BA03408	Medium, Dia 20mm, Hole Dia, 0.8mm
BA03412	Large, Dia 30mm, Hole Dia, 1.2mm



Wire Tightener

Product Code	Specification
BA035260	Length 260mm



Wire Twister/Cutter

Product Code	Specification
BA036150	Length 150mm



Kirschner Wire Punch

Product Code	Specification
BA03701	For 0.6-1.6mm, Length 135mm
BA03702	For 1.5-2.5mm, Length 135mm



K Wire Bender

Product Code	Specification
BA038140	For 0.6-2.5mm, Length 140mm



Bone Curette

Product Code	Specification
BA03902	2mm
BA03903	3mm
BA03904	4mm



Bone Hammer

Product Code	Specification
BA0400105	50g
BA0400214	140g
BA0400327	270g
BA0400435	350g



BW 1.1/1.5/2.0 Reconstruction Plates Bender

Product Code	Specification
BA04101	Length 140mm



2.4/2.7 Reconstruction Plates Bender

Product Code	Specification
BA04201140	Length 140mm



BW 3.5 Reconstruction Plates Bender

Product Code	Specification
BA04302	Length 140mm



Plates Bender

Product Code	Specification
BA04402175	AO 2.7/3.5, Length 175mm



3.5mm Reconstruction Plates Bender

Product Code	Specification
BA04501	Length 280mm

» ORTHOPEDIC INSTRUMENTS




AO Reconstruction Plate Plier

Product Code

BA04601

Specification

Used for 4-9mm width reconstruction plates, including ALPS and MCP, in plane and out of plane bending



Veterinary Plate Model

Product Code

BA04701
BA04702
BA04703
BA04704
BA04705

Specification

Width 4.8mm, apply to 1.5/2.0mm LC-DCP&LCP
Width 6mm, apply to 2.4mm LC-DCP&LCP
Width 8mm, apply to 2.7mm LC-DCP&LCP
Width 10mm, apply to 3.5mm LC-DCP&LCP
Width 12mm, apply to 3.5mm LC-DCP&LCP



Ruskin Ronger

Product Code

BA0480102
BA0480203

Specification

Curved 2mm, Length 160mm
Curved 3mm, Length 160mm



Horsley Bone Cutter

Product Code

BA0490135

Specification

Curved Straight, Length 150mm



Cottle Bone Cutter

Product Code

BA0500125

Specification

Straight 25mm, Length 190mm



Bone Cutter

Product Code

BA0510105

Specification

Length 150mm



AO Type Periosteal Elevator

Product Code

BA0530106
BA0530210

Specification

Flat 6mm, Length 180mm
Curve 10mm, Length 180mm

Double Head Periosteal Elevator

Product Code

BA0520103
BA0520105

Specification

3mm
5mm

» ORTHOPEDIC INSTRUMENTS



Bone Chisel

Product Code

BA05401
BA05402
BA05403
BA05404
BA05405
BA05406
BA05407
BA05408

Specification

3mm
4mm
5mm
6mm
8mm
10mm
12mm
15mm



Bone Rasp

Product Code

BA05501
BA05502

Specification

Curve 5mm Width
Straight 5mm Width



Bone Saw

Product Code

BA05601

Specification

0.4mm Cut, Length 140mm



Trochlea Rasp

Product Code

BA05701
BA05702

Specification

Mini Dia 2.2mm
Small Dia 3.5mm



Double Head Bone Rasp

Product Code

BA05801



Mini Hoffman Retractor

Product Code

BA0590106
BA0590208
BA0590312

Specification

6mm Wide, Length 165mm
8mm Wide, Length 165mm
14mm Wide, Length 165mm



Femoral Bone Retractor

Product Code

BA06001
BA06002
BA06003

Specification

Small 125mmx8mmx12mm
Medium 175mmx13mmx18mm
Large 255mmx25mmx36mm



Skin Hook

Product Code

BA06101



Double Head Hook

Product Code	Specification
BA06201	Sharp Teeth
BA06202	Blunt Teeth



Gelpi Self Retaining Retractor

Product Code	Specification
BA06301	Length 100mm
BA06302	Length 140mm
BA06303	Length 160mm



Gelpi Self Retaining Retractor

Product Code	Specification
BA06401	Length 140mm
BA06402	Length 160mm



Self-Retaining Retractor

Product Code	Specification
BA06501140	Length 140mm
BA06502160	Length 160mm



Self-Retaining Retractor(Point)

Product Code	Specification
BA06601140	Length 140mm
BA06602160	Length 160mm



Orthopedic Measuring Ruler

Product Code	Specification
BA06701	Measuring Length: 150mm
	Measuring Screw Length: 0-60mm
	Measuring Diameter: 0.8/1.0/1.2/1.5/2.0/2.4/2.5/2.7/3.0/3.2/3.5/4.0/4.5/5.0/5.5



Joint Motion Fixator

Product Code	Specification
6100-0101	Ø2.0x70mm 4Pcs
Matching Instruments	
	3.0mmL Type hex Wrench 2.0mm Screw Driver



Straight Fixator

Product Code	Specification
6100-0102	Ø2.0x70mm 4Pcs
Matching Instruments	
	3.0mmL Type hex Wrench 2.0mm Screw Driver



Universal Compression Extension Fixator

Product Code	Specification
6100-0103	Ø2.5x50mm 4Pcs
Matching Instruments	
	3.0mmL Type hex Wrench 2.0mm Screw Driver



Universal Fixator

Product Code	Specification
6100-0104	Ø2.0x50mm 4Pcs
Matching Instruments	
	3.0mmL Type hex Wrench 2.0mm Screw Driver



Combined Fixator

Product Code	Specification
6100-0105	Ø2.0x70mm 4Pcs
Matching Instruments	
	3.0mmL Type hex Wrench 2.0mm Screw Driver



Straight Track Fixator

Product Code	Specification
6100-0106	Ø2.0x70mm 4Pcs
Matching Instruments	
	3.0mmL Type hex Wrench 2.0mm Screw Driver



Horizontal Axis Fixator

Product Code	Specification
6100-0107	Ø2.0x70mm 4Pcs
Matching Instruments	
	3.0mmL Type hex Wrench 2.0mm Screw Driver



Track Vertical Axis Fixator

Product Code	Specification
6100-0108	Ø2.0x70mm 4Pcs
Matching Instruments	
	3.0mmL Type hex Wrench 2.0mm Screw Driver

EXTERNAL FIXATION SYSTEM

Pin to Rod Coupling CB01



Product Code	Specification	Material
6100-0109	Ø3/2	Aluminum Alloy

Note: Can be combined as required

Rod to Rod Coupling CB02



Product Code	Specification	Material
6100-0110	Ø3/3	Aluminum Alloy

Note: Can be combined as required

Rod to Rod Coupling CB03



Product Code	Specification	Material
6100-0111	Ø6/6	Stainless Steel

Note: Can be combined as required

Pin to Rod Coupling CB04




Product Code	Specification	Material
6100-0112	Ø2.5-3/6	Stainless Steel

Note: Can be combined as required

EXTERNAL FIXATION SYSTEM

Schanz Screw SS/Ti



Product Code	DxL(mm)	HA	Diameter(mm)	Length(mm)	Material
6100-0127	Ø2.0x50	HA	Ø2.0	16	Stainless Steel Titanium Alloy
6100-0128	Ø2.0x60			16	
6100-0129	Ø2.0x70			16	
6100-0130	Ø2.5x60	HA	Ø2.5	16	
6100-0131	Ø2.5x70			20	
6100-0132	Ø2.5x80			20	
6100-0133	Ø2.5x90			25	
6100-0134	Ø2.5x100			30	
6100-0135	Ø3.0x60	HA	Ø3.0	20	
6100-0136	Ø3.0x70			20	
6100-0137	Ø3.0x80			25	
6100-0138	Ø3.0x90			25	
6100-0139	Ø3.0x100			30	
6100-0140	Ø3.0x110			35	
6100-0141	Ø3.0x120			35	
6100-0142	Ø3.5x60	HA	Ø3.5	20	
6100-0143	Ø3.5x70			20	
6100-0144	Ø3.5x80			25	
6100-0145	Ø3.5x90			25	
6100-0146	Ø3.5x100			30	
6100-0147	Ø3.5x110			35	
6100-0148	Ø3.5x120			40	

Connecting Rod CD001



Product Code	Diameter(mm)	Length(mm)	Material
6100-0113	3.0	60	Aluminum Alloy
6100-0114	3.0	80	
6100-0115	3.0	100	
6100-0116	6.0	100	
6100-0117	6.0	150	
6100-0118	6.0	200	Carbon Fibre
6100-0119	6.0	250	
6100-0120	3.0	60	
6100-0121	3.0	80	
6100-0122	3.0	100	
6100-01123	6.0	100	
6100-01124	6.0	150	
6100-01125	6.0	200	
6100-01126	6.0	250	

Multi-Functional Power Tools



No.	Product Code	Description	Parameter
1	DA001	Single trigger modular handpiece x1	
2	DB001	Oscillating saw adaptor x1	
3	DB002	AO quick coupling adaptor X1	
4	DB003	Drill bit adaptor x1	
5	DB004	Neuro burs adaptor x1	Short
6	DB005	Neuro burs adaptor x1	Long
7	DB006	K-wire adaptor x1	0.6-2.0mm
8	DB007	TPLO saw adaptor x1	
9	DC001	Battery compartment x1	
10	DC002	Battery x1	
11	DC003	Aseptic transfer kit x1	
12	DC004	Charge x1	
13	DD001	Sterilization box x1	



Cannulated Bone Drill M-06

Quick Coupling Chuck
 Maximum speed: 1200 rpm
 Sterilizable, Positive and negative rotation
 Maximum cannulated diameter: 4.5mm
 Gripping drill bit 0.6-8mm diameter
 Battery voltage: 14.4V
 Charger voltage: 110V/220V



Oscillating Saw M-07

Reciprocating Frequency: 14000 Times/min, 5 degree
 Battery voltage: 14.4V
 Charger voltage: 110V/220V
 Sterilizable
 Blade size: 1.0x24x90mm, 1.0x24x110mm,
 1.2x24x90mm, 1.2x24x110mm,
 1.2x20x120mm



Small Bone Drill M-11

Maximum speed: 1200rpm
 Sterilizable, Positive and negative rotation
 Gripping drilling bit: 0.4-4mm diameter
 Battery voltage: 7.2V
 Charger voltage: 110V/220V



Small Oscillating Saw M-12

Oscillating Frequency: 14000 Times/min, 5 degree
 Battery voltage: 7.2V
 Charger voltage: 110V/220V
 Saw Blade: 0.6x5x70mmx1, 0.6x8x70mmx2,
 0.6x12x70mmx2
 1.2x20x120mm
 Sterilizable



Small Cannulated Bone Drill M-20

Maximum speed: 1200rpm
 Sterilizable, Positive & negative rotation
 Gripping drilling bit: 0.4-4mm diameter
 Battery voltage: 7.2V
 Charger voltage: 110V/220V



Sterilization Box M-21

POWER TOOL SYSTEM

Mini Multi-Functional Drill M-16

- 1- Hand Piece Voltage 80V
- 2- Drill Chuck Maximum speed: 1100rpm, stepless speed regulation, rotating in positive and negative direction for selection, stable operation
- 3- Oscillating Saw Chuck, Reciprocating Frequency: 16000 Times/min;
- 4- AO Chuck Connect all AO quick coupling drill bit, Maximum speed: 1100rpm
- 5- K-wire Chuck Maximum speed: 1100rpm



Mini AO Chuck
Maximum Speed: 1100rpm

Mini Drill Chuck
Maximum Speed: 1100rpm

K-wire Chuck
Maximum Speed: 1100rpm

Mini Oscillating Saw Chuck
Reciprocating Frequency: 16000 Times/Min

Handpiece

TPLO Saw M-17




No.	Product Code	Description	Parameter
1	DA004	TPLO Saw handpiece X1	0-16000cpm
2	DA002	Hexagon Spanner X1	
3	DC001	Battery X2	
4	DC003	Sterilizing Channel X2	
5	DC004	Charger X1	
6	DC005	Aluminum Box X1	
7	DE001	Curved Saw Blade X2	
8	DF001	Sterilization Box X1	

Curved Saw Blade




No.	Product Code	Description
1	M-17-1	R:10mm, 24mmx0.6mm
2	M-17-2	R:12mm, 30mmx0.6mm
3	M-17-3	R:15mm, 30mmx0.6mm
4	M-17-4	R:18mm, 35mmx0.6mm
5	M-17-5	R:20mm, 35mmx0.6mm
6	M-17-6	R:24mm, 42mmx0.6mm
7	M-17-7	R:27mm, 50mmx0.6mm
8	M-17-8	R:30mm, 50mmx0.6mm



Sterilization Box

Product Code DC006
Specification 250x220x80mm



Sterilization Box for TPLO Saw


Product Code DC007
Specification 300x245x80mm

GENERAL SURGERY INSTRUMENTS



Mayo Safety Pin

Product Code EA001
Specification 140mm



Halsted-Mosquito

Product Code EA00101
EB00102
Specification 125mm Straight
125mm Curved



Crile Artery Forceps

Product Code EB00201
EB00202
EB00203
EB00204
Specification 140mm Straight
140mm Curved
160mm Straight
160mm Curved




Allis Tissue Forceps

Product Code EB00301
EB00302
Specification 125mm/4x5 teeth
150mm/4x5 teeth



Dressing Blunt

Product Code EC00101
EC00102
EC00103
EC00104
Specification 145mm Straight
145mm Curved
165mm Straight
165mm Curved



Dressing Sharp

Product Code EC00201
EC00202
EC00203
EC00204
EC00205
EC00206
Specification 115mm Straight
115mm Curved
145mm Straight
145mm Curved
165mm Straight
165mm Curved



Thread Cutter

Product Code	Specification
EC003	/



Pin Holder With Cutter

Product Code	Specification
EC00401	140mm
EC00402	160mm



Dissecting Forceps

Product Code	Specification
ED001	150mm



Adson Dissecting Forceps

Product Code	Specification
ED002	With Jaw



Adson Dissecting Forceps

Product Code	Specification
ED003	With Hook



Adson Dissecting Forceps

Product Code	Specification
ED004	Brown



Standard Scalpel Handles

Product Code	Specification
ED00501	#3
ED00502	#4



Backhuas Towel Forces

Product Code	Specification
EE00101	100mm
EE00102	125mm



Backhuas Towel Forces

Product Code	Specification
EE002	90mm



Spay Ligation Fixture Clamp

Product Code	Specification
EE00301	140mm
EE00302	160mm




Canine/Felines Snook Hook

Product Code	Specification
EF00101	150mm felines
EF00102	160mm canines



Tongue Forceps

Product Code	Specification
EG001	170mm



Balfour Abdominal Retractor

Product Code	Specification
EH001	Max 180mm

Notes
can be used as the rib retractor, open abdominal surgery hook.



Tail Cutter

Product Code	Specification
EI001	Length 140mm




Mouth Opener

Product Code	Specification
EJ001	100mm, range 0-25mm


» GENERAL SURGERY INSTRUMENTS

Veterinary Soft Tissue Surgical Orthopedic Instrument Set E001



No.	Product Code	Item	Specification
1	E001-01	Mayo safety pin X1	
2	E001-02	Crile Artery Forceps(Straight) X3	160mm
3	E001-03	Crile Artery Forceps(Curved) X3	160mm
4	E001-04	Allis Tissue Forceps X2	
5	E001-05	Backhuas Towel Forces X4	125mm
6	E001-06	Dressing Blunil X1	
7	E001-07	Metzenbaum surper Dressing Blunt (Straight) X1	140mm
8	E001-08	Metzenbaum surper Dressing Blunt (Curved) X1	140mm
9	E001-09	Pin Holder With Cutter X1	160mm
10	E001-10	Adson Dissecting Forceps(with jaw) X1	
11	E001-11	Adson Dissecting Forceps(with hook) X1	
12	E001-12	Standard scalpel handles X1	#3/4
13	E001-13	Sterilization case(selective)	

Spay Pack Veterinary Instrument Set E002



No.	Product Code	Item	Specification
1	E002-01	Mayo safety pin X1	
2	E002-02	Mayo Dressing Blunt X1	150mm
3	E002-03	Dressing Blunt (Curved) X1	160mm
4	E002-04	Crile Artery Forceps(Curved) X2	140mm
5	E002-05	Crile Artery Forceps(Straight) X2	160mm
6	E002-06	Crile Artery Forceps(Straight) X3	125mm
7	E002-07	Holsted-mosquito(Curved)x2	125mm
8	E002-08	Backhuas Towel Forces X4	90mm
9	E002-09	Allis Tissue Forceps X2	140mm
10	E002-10	Dissecting Forceps X1	
11	E002-11	Pin Holder With Cutter X1	160mm
12	E002-12	Adson Dissecting Forceps X1	
13	E002-13	Standard scalpel handles X1	#3/4
14	E002-15	Sterilization box	

» DENTAL INSTRUMENTS



Sickle Scaler

Product Code
FA001-01

Specification
Length 170mm



Very Fine Periosteal Elevator

Product Code
FA001-02

Specification
Length 170mm, 2/4mm



Subgingival Curette

Product Code
FA001-03

Specification
Length 170mm



Jaquette Subgingival Scaler

Product Code
FA001-04

Specification
Length 170mm



Sickle Scaler

Product Code
FA001-05

Specification
Length 170mm



Tartar Removing Forceps

Product Code
FA001-06

Specification
Length 125mm



Extraction Forceps

Product Code
FA001-07

Specification
Length 125mm



Fragment Extraction Forceps

Product Code
FA001-08

Specification
Length 140mm



Periosteal Elevator

Product Code
FA001-09

Specification
Length 150mm



Straight Root Tip Pick

Product Code
FA001-10

Specification
Length 150mm



Winged Elevator

Product Code
FA001-11
FA001-12
FA001-13
FA001-14

Specification
Head width 1mm, Length 150mm
Head width 2mm, Length 150mm
Head width 3mm, Length 150mm
Head width 4mm, Length 150mm





Minnesota Retractor

Product Code
FA001-15

Specification
155mm

Dental Instrument Kit FA001



No.	Product Code	Description	Specification
1	FA001-01	Sickle Scaler	170mm
2	FA001-02	Very Fine Periosteal Elevator(sharp)	170mm, 2/4mm
3	FA001-03	Subgingivar Curette	170mm
4	FA001-04	Joquette Supragingival Scaler (with measure)	170mm
5	FA001-05	Sickle Scaler	170mm
6	FA001-06	Tartar Removing Forceps	125mm
7	FA001-07	Extraction Forceps	125mm
8	FA001-08	Fragment Extraction Forceps	140mm
9	FA001-09	Periosteal Elevator	150mm
10	FA001-10	Straight Rool Tip Pick	150mm
11	FA001-11	Winged Elevator	1mm, 150mm
12	FA001-12	Winged Elevator	2mm, 150mm
13	FA001-13	Winged Elevator	3mm, 150mm
14	FA001-14	Winged Elevator	4mm, 150mm
15	FA001-15	Minnesota Retractor	155mm



Orthopromed Inc.

Address: 964 E. Badillo Street #2140
Covina CA 91724 **UNITED STATES**

Phone: 302-613-4333

Web: www.orthopromed.com

Email: info@orthopromed.com

After Sales: support@orthopromed.com

